

Datasheet for ABIN3115509

## FAM75A6 Protein (AA 1-1343) (Strep Tag)



[Go to Product page](#)

### Overview

|                               |                                                  |
|-------------------------------|--------------------------------------------------|
| Quantity:                     | 250 µg                                           |
| Target:                       | FAM75A6 (SPATA31A6)                              |
| Protein Characteristics:      | AA 1-1343                                        |
| Origin:                       | Human                                            |
| Source:                       | Cell-free protein synthesis (CFPS)               |
| Protein Type:                 | Recombinant                                      |
| Purification tag / Conjugate: | This FAM75A6 protein is labelled with Strep Tag. |
| Application:                  | ELISA, SDS-PAGE (SDS), Western Blotting (WB)     |

### Product Details

|           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Brand:    | AliCE®                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Sequence: | <p>MENLPFPLKL LSASSLNAPS STPWVLDIFL TLVFALGFFF LLLPYLSYFH CDDPPSPSPG</p> <p>KRKCPVGRRR RPRGRMKNHS LRAGRECPRG LEETSDLLSQ LQSLLGPHLD KGDFGQLSGP</p> <p>DPPGEVGERA PDGASQSSHE PMEDAAPILS PLASPDPAK HPQDLASTPS PGPMTTSVSS</p> <p>LSASQPPEPS LPLEHPSPEP PALFPHPPHT PDPLACSPPP PKGFTAPPLR DSTLITPSHC</p> <p>DSVALPLGTV PQSLSPHEDL VASVPAISGL GGSNSHVSAS SRWQETARTS CAFNSSVQQD</p> <p>PLSRHPPETC QMEAGSLFLL SSDGQNVVGI QVTETAKVNI WEEKENVGSF TNQMTPEKHL</p> <p>NSLGNLAKSL DAEQDTTNPK PFWNMGENSK QLPGPQKCSD PRLLQESFWK NYSQFLWGLP</p> <p>SLHSESLVAN AWTDRSYTL QSPFLFNEM SNVCPIQRET TMSPLLFQAQ PLSHRQPFIS</p> <p>STPQFLPTPM AQAEQAHLQ SSFPVLSPAF PSLIKNTGVA CPASQNKVQA LSLPETQHPE</p> <p>WPLLRLKQLEG RLALPSRVQK SQDVFSVSTP NLPQESLTSI LPENFPVSPE LRRQLEQHIK</p> <p>KWIIQHWGNL GRIQESLDLM QLRDESPGTS QAKGKPSPWQ SSTSTGESSK EAQKVQFQLE</p> |

RDLCPHLGQI LGETPQNLSR DMKSFPRKVL GVTSEESERN LRKPLRSDSG SDLLRCTERT  
HIENILKAHM GRNLGQTNEG LIPVRVRRSW LAVNQALPVS NTHVKTSNLA APKSGKACVN  
TAQVLSFLEP CTQQGLGAHI VRFWAKHRWG LPLRVLKPIQ CFKLEKVSSL SLTQLAGPSS  
ATCESGAGSE VEVDMLRKP PMASLRKQVL TKASDHMPES LLASSPAWKQ FQRAPRGIPS  
WNDHGPKLPP PAGQEGRWPS KPLTYSLTGS TQQSRSLGAQ SSKAGETREA VPQCRVPLET  
CMLANLQATS EDVHGFEAPG TSKSSLHPRV SVSQDPRKLC LMEEVVSEFE PGMATKSETQ  
PQVCAAVVLL PDGQASVVPV ASENLVSQVP QGHLQSMPTG NMRASQELHD LMAARRSKLV  
QEEPRNPNCQ GSCKSQRPMF PPIHKSEKSR KPNLEKHEER LEGLRTPQLT PVRKTEDTHQ  
DEGVQLPSK KQPPSVSHFG ENIKQFFQWI FSKKKSKPAP VTAESQKTVK NRSCVYSSSA  
EAQGLMTAVG QMLDKKMSLC HAHHASKVNQ HKQKFQAPVC GFPCNHRHLF YSEHGRILSY  
AASSQQATLK SQGCPNRDRQ IRNQQLPSV RCNNEQWGLR HPQILHPKKA VSPVSPQHW  
PKTSGASSHH HHCPRHCLLW EGI

**Sequence without tag. The proposed Strep-Tag is based on experience s with the expression system, a different complexity of the protein could make another tag necessary. In case you have a special request, please contact us.**

---

### Characteristics:

#### Key Benefits:

- Made in Germany - from design to production - by highly experienced protein experts.
- Protein expressed with ALiCE® and purified in one-step affinity chromatography
- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab try to ensure that you receive soluble protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

#### Expression System:

- ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from *Nicotiana tabacum* c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications.
- During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to

## Product Details

produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!

### Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured against its specific reference buffer.
- We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein.

|               |                                                                                                              |
|---------------|--------------------------------------------------------------------------------------------------------------|
| Purification: | One-step Strep-tag purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®). |
| Purity:       | > 70-80 % as determined by SDS PAGE, Western Blot and analytical SEC (HPLC).                                 |
| Grade:        | custom-made                                                                                                  |

## Target Details

|                   |                                                                                                                        |
|-------------------|------------------------------------------------------------------------------------------------------------------------|
| Target:           | FAM75A6 (SPATA31A6)                                                                                                    |
| Alternative Name: | SPATA31A6                                                                                                              |
| Background:       | Spermatogenesis-associated protein 31A6 (Protein FAM75A6),FUNCTION: May play a role in spermatogenesis. {ECO:0000250}. |
| Molecular Weight: | 147.8 kDa                                                                                                              |
| UniProt:          | <a href="#">Q5VVP1</a>                                                                                                 |

## Application Details

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Application Notes: | In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.                                                                                                                                                                                                                                                                                                                                                                                        |
| Comment:           | <p>ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from <i>Nicotiana tabacum</i> c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications.</p> <p>During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional</p> |

## Application Details

components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!

Restrictions: For Research Use only

## Handling

Format: Liquid

Buffer: The buffer composition is at the discretion of the manufacturer.  
Standard Storage Buffer: PBS pH 7.4, 10 % Glycerol **Might differ depending on protein.**

Handling Advice: Avoid repeated freeze-thaw cycles.

Storage: -80 °C

Storage Comment: Store at -80°C.

Expiry Date: 12 months