antibodies

Datasheet for ABIN3117454 SLC6A18 Protein (AA 1-628) (Strep Tag)





Overview

| Quantity: | 1 mg |
|-------------------------------|--|
| Target: | SLC6A18 |
| Protein Characteristics: | AA 1-628 |
| Origin: | Human |
| Source: | Tobacco (Nicotiana tabacum) |
| Protein Type: | Recombinant |
| Purification tag / Conjugate: | This SLC6A18 protein is labelled with Strep Tag. |
| Application: | ELISA, SDS-PAGE (SDS), Western Blotting (WB) |

Product Details

| Sequence: | MAHAPEPDPA ACDLGDERPK WDNKAQYLLS CTGFAVGLGN IWRFPYLCQT YGGGAFLIPY |
|-----------|---|
| | VIALVFEGIP IFHVELAIGQ RLRKGSVGVW TAISPYLSGV GLGCVTLSFL ISLYYNTIVA |
| | WVLWYLLNSF QHPLPWSSCP PDLNRTGFVE ECQGSSAVSY FWYRQTLNIT ADINDSGSIQ |
| | WWLLICLAAS WAVVYMCVIR GIETTGKVIY FTALFPYLVL TIFLIRGLTL PGATKGLIYL |
| | FTPNMHILQN PRVWLDAATQ IFFSLSLAFG GHIAFASYNS PRNDCQKDAV VIALVNRMTS |
| | LYASIAVFSV LGFKATNDYE HCLDRNILSL INDFDFPEQS ISRDDYPAVL MHLNATWPKR |
| | VAQLPLKACL LEDFLDKSAS GPGLAFVVFT ETDLHMPGAP VWAMLFFGML FTLGLSTMFG |
| | TVEAVITPLL DVGVLPRWVP KEALTGLVCL VCFLSATCFT LQSGNYWLEI FDNFAASPNL |
| | LMLAFLEVVG VVYVYGMKRF CDDIAWMTGR RPSPYWRLTW RVVSPLLLTI FVAYIILLFW |
| | KPLRYKAWNP KYELFPSRQE KLYPGWARAA CVLLSLLPVL WVPVAALAQL LTRRRRTWRD |
| | RDARPDTDMR PDTDTRPDTD MRPDTDMR |
| | Sequence without tag. The proposed Strep-Tag is based on experience s with the expression |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3117454 | 04/16/2024 | Copyright antibodies-online. All rights reserved.

| | system, a different complexity of the protein could make another tag necessary. In case you |
|------------------|---|
| | have a special request, please contact us. |
| Characteristics: | Key Benefits: |
| | Made in Germany - from design to production - by highly experienced protein experts. Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure correct folding and modification. These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed). State-of-the-art algorithm used for plasmid design (Gene synthesis). |
| | This protein is a made-to-order protein and will be made for the first time for your order. Our |
| | experts in the lab will ensure that you receive a correctly folded protein. |
| | The big advantage of ordering our made-to-order proteins in comparison to ordering custom |
| | made proteins from other companies is that there is no financial obligation in case the protein |
| | cannot be expressed or purified. |
| | Expression System: |
| | ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post- translational modifications. |
| | During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein! |
| | Concentration: |
| | The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer. We use the Eveness's Protection tool to determine the characteristic coefficient of each protein. |
| | • We use the Expasy's ProtParam tool to determine the absorption coefficient of each protein. |
| Purification: | Two step purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®): |
| | 1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3117454 | 04/16/2024 | Copyright antibodies-online. All rights reserved. capture material. Eluate fractions are analyzed by SDS-PAGE.2. Protein containing fractions of the best purification are subjected to second purification step
through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and
Western blot.Purity:>80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.Endotoxin Level:Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)</td>Grade:Crystallography grade

Target Details

| Target: | SLC6A18 |
|-------------------|--|
| Alternative Name: | SLC6A18 (SLC6A18 Products) |
| Background: | Inactive sodium-dependent neutral amino acid transporter B(0)AT3 (Sodium- and chloride- dependent transporter XTRP2) (Solute carrier family 6 member 18) (System B(0) neutral amino acid transporter AT3),FUNCTION: Does not show neutral amino acid transporter activity. {ECO:0000269 PubMed:26240152}. |
| Molecular Weight: | 70.9 kDa |
| UniProt: | Q96N87 |

Application Details

| Application Notes: | In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though. |
|--------------------|--|
| Comment: | ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications. During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein! |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN3117454 | 04/16/2024 | Copyright antibodies-online. All rights reserved.

Application Details

Restrictions:

For Research Use only

Handling

| Format: | Liquid |
|------------------|--|
| Buffer: | The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us. |
| Handling Advice: | Avoid repeated freeze-thaw cycles. |
| Storage: | -80 °C |
| Storage Comment: | Store at -80°C. |
| Expiry Date: | Unlimited (if stored properly) |

Images




Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process