

Datasheet for ABIN3117494

**SORCS2 Protein (AA 51-1159) (rho-1D4 tag)**[Go to Product page](#)**1** Image

## Overview

Quantity:	1 mg
Target:	SORCS2
Protein Characteristics:	AA 51-1159
Origin:	Human
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This SORCS2 protein is labelled with rho-1D4 tag.
Application:	ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS)

## Product Details

Sequence:	RSPEPGRGLGP HAQLTRVPRS PPAGRAEPGG GEDRQARGTE PGAPGPSPGP APGPGEDGAP AAGYRRWERA APLAGVASRA QVSLISTSFV LKGDATHNQA MVHWTGENSS VILILTKYYH ADMGKVLESS LWRSSDFGTS YTKLTLQPGV TTVIDNFYIC PTNKRKVILV SSSLSDRDQS LFLSADEGAT FQKQPIPFV ETLIFHPKEE DKVLAYTKES KLYVSSDLGK KWTLLQERVT KDHVFWVSVG VDADPDLVHV EAQDLGGDFR YVTCAIHNCES EKMLTAPFAG PIDHGSLTVQ DDYIFFKATS ANQTKYYVSY RRNEFVLMKL PKYALPKDLQ IISTDESQVF VAVQEWWQMD TYNLYQSDPR GVRVALVLQD VRSSRQAEES VLIDILEVRG VKGVFLANQK IDGKVMTLIT YNKGRDWDYL RPPSMDMNGK PTNCKPPDCH LHLHLRWADN PYVSGTVHTK DTAPGLIMGA GNLGSQVLEY KEEMYITSDC GHTWRQVFEE EHHILYLDHG GVIVAIKDTS IPLKILKFSV DEGLTWSTHN FTSTSVFVDG LLSEPGDETL VMTVFGHISF RSDWELVKVD FRPSFSRQCG EEDYSSWELS NLQGDRICMG QQRSFRKRKS TSWCIKGRSF TSALTSRVCE CRDSDFLCDY GFERSSSES STNKCSANFW FNPLSPDDC ALGQTYTSSL GYRKVVSNCV EGGVDMQQSQ
-----------	---

VQLQCPLTPP RGLQVSIQGE AVAVRPGEDV LFVVRQEQGD VLTKYQVDL GDGFKAMYVN  
LTLTGEPRIH RYESPGIYRV SVRAENTAGH DEAVLFVQVN SPLQALYLEV VPVIGLNQEV  
NLTA VLLPLN PNLTVFYWWI GHSLQPLLSL DNSVTTRFSD TGDVRVTVQA ACGNSVLQDS  
RVLRLVDQFQ VMPLQFSKEL DAYNPNTPEW REDVGLVVTR LLSKETSVPQ ELLVTVVKPG  
LPTLADLYVL LPPRPTRKR SLSSDKRLAA IQQVLNAQKI SFLLRGGVRV LVALRDTGTG  
AEQLGGGGGY WAVVVLVFIG LFAAGAFILY KFKRKRPGRV VYAQMHNKE QEMTSPVSHS  
EDVQGAVQGN HSGVLSINS REMHSYLV

**Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.**

---

### Characteristics:

- Made in Germany - from design to production - by highly experienced protein experts.
- Human SORCS2 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receipt of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the ExPASy's protparam tool to determine the absorption coefficient of each protein.

---

### Purification:

Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:

1. Membrane proteins are fractionated by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.
2. The best performing detergent is used for solubilization and the proteins are purified via their

## Product Details

rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.

3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility:	0.22 µm filtered
Endotoxin Level:	Protein is endotoxin-free.
Grade:	Crystallography grade

## Target Details

Target:	SORCS2
Alternative Name:	SORCS2 ( <a href="#">SORCS2 Products</a> )
Molecular Weight:	124.4 kDa Including tag.
UniProt:	<a href="#">Q96PQ0</a>

## Application Details

Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only

## Handling

Format:	Liquid
Buffer:	100 mM NaCl, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C

## Handling

---

Storage Comment: Store at -80°C.

Expiry Date: Unlimited (if stored properly)

## Images

---



**Image 1.** „Crystallography Grade“ protein due to multi-step, protein-specific purification process