antibodies -online.com





SCN11A Protein (AA 1-126) (His tag)



Image



Go to Product page

\sim	
()\/△	rview
\cup	1 410 44

s labelled with His tag.
g (WB), Crystallization (Crys), SDS-PAGE (SDS)
RPFT FDSLAAIEKR ITIQKEKKKS KDKAATEPQP RPQLDLKASR
PLED LDPFYKDHKT FMVLNKKRTI YRFSAKRALF ILGPFNPIRS FMIRIS
. Tag location is at the discretion of the manufacturer. If you have a
e contact us.
from design to production - by highly experienced protein experts. ein (raised in Insect Cells) purified by multi-step, protein-specific process tion grade. rithm used for plasmid design (Gene synthesis).
to order protein and will be made for the first time for your order. Our
ensure that you receive a correctly folded protein.
ordering our made-to-order proteins in comparison to ordering custom
her companies is that there is no financial obligation in case the protein
F

cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

SCN11A

Endotoxin Level:

Protein is endotoxin free.

Grade:

Target:

Crystallography grade

Target Details

Alternative Name:	Scn11a (SCN11A Products)
Background:	This protein mediates the voltage-dependent sodium ion permeability of excitable membranes.
	Assuming opened or closed conformations in response to the voltage difference across the
	membrane, the protein forms a sodium-selective channel through which sodium ions may pass
	in accordance with their electrochemical gradient. It is a tetrodotoxin-resistant sodium channel
	isoform. Also involved, with the contribution of the receptor tyrosine kinase NTRK2, in rapid

Storage Comment:

Expiry Date:

Store at -80°C.

Unlimited (if stored properly)

Target Details	
	BDNF-evoked neuronal depolarization (By similarity). {ECO:0000250}.
Molecular Weight:	15.9 kDa Including tag.
UniProt:	Q9R053
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though.
Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C



Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process