antibodies - online.com







ATP5E Protein (AA 2-52) (His tag)





\sim							
	1//	\Box	$r \setminus$	/ [\bigcirc	1	٨,

1 mg
ATP5E
AA 2-52
Mouse
Escherichia coli (E. coli)
Recombinant
This ATP5E protein is labelled with His tag.
ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS)
VAYWRQAGLS YIRFSQICAK AVRDALKTEF KANAEKTSGS SIKIVKVSKK E
Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a
special request, please contact us.
 Made in Germany - from design to production - by highly experienced protein experts. Mouse Atp5e Protein (raised in E. Coli) purified by multi-step, protein-specific process to ensure crystallization grade.
State-of-the-art algorithm used for plasmid design (Gene synthesis).

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in bacterial culture:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Endotoxin has not been removed. Please contact us if you require endotoxin removal.

Grade:

Crystallography grade

Target Details

Target:	ATP5E
Alternative Name:	Atp5e (ATP5E Products)
Background:	Mitochondrial membrane ATP synthase (F(1)F(0) ATP synthase or Complex V) produces ATP
	from ADP in the presence of a proton gradient across the membrane which is generated by
	electron transport complexes of the respiratory chain. F-type ATPases consist of two structural
	domains, $F(1)$ - containing the extramembraneous catalytic core, and $F(0)$ - containing the
	membrane proton channel, linked together by a central stalk and a peripheral stalk. During
	catalysis, ATP synthesis in the catalytic domain of F(1) is coupled via a rotary mechanism of

Target Details

Taryer Details	
	the central stalk subunits to proton translocation. Part of the complex F(1) domain and of the central stalk which is part of the complex rotary element. Rotation of the central stalk against the surrounding alpha(3)beta(3) subunits leads to hydrolysis of ATP in three separate catalytic sites on the beta subunits (By similarity). {ECO:0000250}.
Molecular Weight:	6.7 kDa Including tag.
UniProt:	P56382
Pathways:	Proton Transport, Ribonucleoside Biosynthetic Process
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though.
Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)



Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process