-online.com antibodies

Datasheet for ABIN3124858 OSCP1 Protein (AA 1-379) (Strep Tag)



0	•	
$() \vee \Box$	r_{1}	
016	rview	

Quantity:1 mgTarget:OSCP1Protein Characteristics:AA 1-379Origin:MouseSource:Tobacco (Nicotiana tabacum)Protein Type:RecombinantProtein tag / Conjugate:This OSCP1 protein is labelled with Strep Tag.Application:SDS-PAGE (SDS), Western Blotting (WB), ELISAProduct DetailsSequence:MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expressior system, a different complexity of the protein could make another tag necessary. In case you		
Protein Characteristics: AA 1-379 Origin: Mouse Source: Tobacco (Nicotiana tabacum) Protein Type: Recombinant Purification tag / Conjugate: This OSCP1 protein is labelled with Strep Tag. Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AikGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Quantity:	1 mg
Origin: Mouse Source: Tobacco (Nicotiana tabacum) Protein Type: Recombinant Purification tag / Conjugate: This OSCP1 protein is labelled with Strep Tag. Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Target:	OSCP1
Source: Tobacco (Nicotiana tabacum) Protein Type: Recombinant Purification tag / Conjugate: This OSCP1 protein is labelled with Strep Tag. Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLUMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL	Protein Characteristics:	AA 1-379
Protein Type: Recombinant Purification tag / Conjugate: This OSCP1 protein is labelled with Strep Tag. Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Origin:	Mouse
Purification tag / Conjugate: This OSCP1 protein is labelled with Strep Tag. Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Source:	Tobacco (Nicotiana tabacum)
Application: SDS-PAGE (SDS), Western Blotting (WB), ELISA Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Protein Type:	Recombinant
Product Details Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Purification tag / Conjugate:	This OSCP1 protein is labelled with Strep Tag.
Sequence: MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Application:	SDS-PAGE (SDS), Western Blotting (WB), ELISA
LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Product Details	
AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression	Sequence:	MSVRTLPLLF LNLGGEMLYV LDQRLRAQNI PGDKARKVLN DIISTMFNRK FMDELFKPQE
VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression		LYSKKALRTV YDRLAHASIM RLNQASMDKL YDLMTMAFKY QVLLCPRPKD VLLVTFNHLD
DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression		AIKGFVQDSP TVIHQVDETF RQLSEVYGKL SEGEFQLIRQ TLLNFFQDLH IRVSTFLKDK
PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression		VQNSNGRFVL PVSGPVPWGT EVPGVIRVFS VKGKEVKKMK FRHGGDYVAA QKEGSFELYG
QLEQNTSKGD DLLAMMDRL Sequence without tag. The proposed Strep-Tag is based on experience s with the expression		DRVLKLGTNM YSASRPVETH MSATSKNAAS RAQENIVPNP LAKEELNFLA RLMGGMEIKK
Sequence without tag. The proposed Strep-Tag is based on experience s with the expression		PSGPEPGFRL NLFTTDEEEE HAALSRPEEL SYEVISIQAT QDQQRNEELA RIMGEFEITE
		QLEQNTSKGD DLLAMMDRL
system, a different complexity of the protein could make another tag necessary. In case you		Sequence without tag. The proposed Strep-Tag is based on experience s with the expression
		system, a different complexity of the protein could make another tag necessary. In case you

Characteristics: Key Benefits:

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/4 | Product datasheet for ABIN3124858 | 05/07/2024 | Copyright antibodies-online. All rights reserved.

have a special request, please contact us.

- · Made in Germany from design to production by highly experienced protein experts.
- Protein expressed with ALICE® and purified by multi-step, protein-specific process to ensure correct folding and modification.
- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

Expression System:

- ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require posttranslational modifications.
- During lysate production, the cell wall and other cellular components that are not required for
 protein production are removed, leaving only the protein production machinery and the
 mitochondria to drive the reaction. During our lysate completion steps, the additional
 components needed for protein production (amino acids, cofactors, etc.) are added to
 produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

Concentration:

- The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in Almost Living Cell-Free Expression System (ALiCE®):

- 1. In a first purification step, the protein is purified from the cleared cell lysate using StrepTag capture material. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/4 | Product datasheet for ABIN3124858 | 05/07/2024 | Copyright antibodies-online. All rights reserved.

Product Details

Purity:	\geq 80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Endotoxin Level:	Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)

Target Details

Target:	OSCP1
Alternative Name:	Oscp1 (OSCP1 Products)
Background:	Protein OSCP1,FUNCTION: May be involved in drug clearance in the placenta. {ECO:0000250}.
Molecular Weight:	43.2 kDa
UniProt:	Q8BHW2

Application Details

Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	 ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications. During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/4 | Product datasheet for ABIN3124858 | 05/07/2024 | Copyright antibodies-online. All rights reserved.

Handling

Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)