antibodies .- online.com





LRFN5 Protein (AA 18-529) (His tag)



Go to Product page

Overview

Quantity:	1 mg
Target:	LRFN5
Protein Characteristics:	AA 18-529
Origin:	Mouse
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This LRFN5 protein is labelled with His tag.
Application:	ELISA, Crystallization (Crys), SDS-PAGE (SDS), Western Blotting (WB)

Product Details

Sequence:

QICPKRCVCQ ILSPNLATLC AKKGLLFVPP NIDRRTVELR LADNFVTNIK RKDFANMTSL VDLTLSRNTI SFITPHAFAD LRNLRALHLN SNRLTKITND MFSGLSNLHH LILNNNQLTL ISSTAFDDVF ALEELDLSYN NLETIPWDAV EKMVSLHTLS LDHNMIDNIP KGTFSHLHKM TRLDVTSNKL QKLPPDPLFQ RAQVLATSGI ISPSTFALSF GGNPLHCNCE LLWLRRLSRE DDLETCASPA LLTGRYFWSI PEEEFLCEPP LITRHTHEMR VLEGQRATLR CKARGDPEPA IHWISPEGKL ISNATRSLVY DNGTLDILIT TVKDTGAFTC IASNPAGEAT QTVDLHIIKL PHLLNSTNHI HEPDPGSSDI STSTKSGSNA SSSNGDTKMS QDKIVVAEAT SSTALLKFNF QRNIPGIRMF QIQYNGTYDD TLVYRMIPPT SKTFLVNNLA SGTMYDLCVL AIYDDGITSL TATRVVGCIQ FTTEQDYVRC HFMQSQFLGG TM

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Lrfn5 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Protein is endotoxin free.

Grade:

Crystallography grade

Target Details

Expiry Date:

LRFN5
Lrfn5 (LRFN5 Products)
Cell adhesion molecule that mediates homophilic cell-cell adhesion in a Ca(2+)-independent manner. Promotes neurite outgrowth in hippocampal neurons (By similarity). {ECO:0000250}.
58.0 kDa Including tag.
Q8BXA0
In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though.
Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
For Research Use only
Liquid
100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Avoid repeated freeze-thaw cycles.
-80 °C
Store at -80°C.

Unlimited (if stored properly)