# antibodies .- online.com





# TRAPPC4 Protein (AA 1-219) (His tag)



## **Image**



( )	1 /	$\sim$	KI /	110	Νę
	1//	$\vdash$	I \/	1 ←	٠// ٢

Overview		
Quantity:	1 mg	
Target:	TRAPPC4	
Protein Characteristics:	AA 1-219	
Origin:	Mouse	
Source:	Escherichia coli (E. coli)	
Protein Type:	Recombinant	
Purification tag / Conjugate:	This TRAPPC4 protein is labelled with His tag.	
Application:	SDS-PAGE (SDS), ELISA, Western Blotting (WB), Crystallization (Crys)	
Product Details		
Sequence:	MAIFSVYVVN KAGGLIYQWD SYSPRAEAEK TFSYPLDLLL KLHDERVLVA FGQRDGIRVG	
	HAVLAINGMD VNGKYTADGK EVLEYLGNPA NYPVSIRFGR PRLTSNEKLM LASMFHSLFA	
	IGSQLSPEQG SSGIEMLETD TFKLHCFQTL TGIKFVVLAD PRQAGIDSLL RKIYEIYSDF	
	ALKNPFYSLE MPIRCELFDQ NLKLALEVAE KAGTFGPGS	
	Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a	
	special request, please contact us.	

#### Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Trappc4 Protein (raised in E. Coli) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

#### Purification:

Two step purification of proteins expressed in bacterial culture:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Endotoxin has not been removed. Please contact us if you require endotoxin removal.

Grade:

Crystallography grade

#### Target Details

Target:	TRAPPC4	
Alternative Name:	Trappc4 (TRAPPC4 Products)	
Background:	May play a role in vesicular transport from endoplasmic reticulum to Golgi. May play a role in dendrite postsynaptic membrane trafficking.	
Molecular Weight:	25.3 kDa Including tag.	

## **Target Details**

UniProt:	Q9ES56
----------	--------

## **Application Details**

Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies		
	as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee		
	though.		
Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the		
	recombinant protein with the default tag will be insoluble our protein lab may suggest a higher		
	molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible		


options with you in detail to assure that you receive your protein of interest.

Restrictions: For Research Use only

# Handling

Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)

### **Images**



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process