

Datasheet for ABIN3131044

Heparanase 2 Protein (HPSE2) (AA 39-592) (His tag)



Overview

Quantity:	1 mg
Target:	Heparanase 2 (HPSE2)
Protein Characteristics:	AA 39-592
Origin:	Mouse
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This Heparanase 2 protein is labelled with His tag.
Application:	ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS)

Product Details

Sequence:

FHAGNRRPLP VDRATGLKEK TLILLDVSTK NPVRTVNENF LSLQLDPSII HDGWLDFLSS

KRLVTLARGL SPAFLRFGGK RTDFLQFQNL RNPAKSRGGP GPDYYLKNYE DDIVRSDVAL

DKQKGCKIAQ HPDVMLELQR EKASQMHLVL LKEQYSNTYS NLILTARSLD KLYNFADCSG

LHLIFALNAL RRNPNNSWNS SSALSLLKYS ASKKYNISWE LGNEPNNYRS IHGRAVNGSQ

LGKDYIQLKS LLQPIRVYSR ASLYGPNIGR PRKNVIALLD GFMKVAGSTV DAVTWQHCYI

DGRVVKVMDF LKTRLLDTLS DQIRKIQKVV NTYTPGKKIW LEGVVTTSAG GTNNLSDSYA

AGFLWLNTLG MLANQGIDVV IRHSFFDHGY NHLVDQNFNP LPDYWLSLLY KRLIGPKVLA

VHVAGLQRKP RPGRVIRDKL RIYAHCTNHH NHNYVRGSIT LFIINLHRSR KKIKLAGTLR

DKLVHQYLLQ PYGQEGLKSK SVQLNGQPLV MVDDGTLPEL KPRPLRAGRT LVIPPVTMGF

YVVKNVNALA CRYR

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Hpse2 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- 2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Protein is endotoxin free.

Grade:

Crystallography grade

Target Details

Expiry Date:

Heparanase 2 (HPSE2)
Hpse2 (HPSE2 Products)
Binds heparin and heparan sulfate with high affinity, but lacks heparanase activity. Inhibits
HPSE, possibly by competing for its substrates (in vitro) (By similarity). {ECO:0000250}.
63.6 kDa Including tag.
B2RY83
Glycosaminoglycan Metabolic Process
In addition to the applications listed above we expect the protein to work for functional studies
as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee
though.
Protein has not been tested for activity yet. In cases in which it is highly likely that the
recombinant protein with the default tag will be insoluble our protein lab may suggest a higher
molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible
options with you in detail to assure that you receive your protein of interest.
For Research Use only
Liquid
100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Avoid repeated freeze-thaw cycles.
-80 °C
Store at -80°C.

Unlimited (if stored properly)