antibodies .- online.com





SERPINA3K Protein (AA 22-418) (His tag)



Image



Overview

Quantity:	1 mg
Target:	SERPINA3K
Protein Characteristics:	AA 22-418
Origin:	Mouse
Source:	Escherichia coli (E. coli)
Protein Type:	Recombinant
Purification tag / Conjugate:	This SERPINA3K protein is labelled with His tag.
Application:	ELISA, Crystallization (Crys), SDS-PAGE (SDS), Western Blotting (WB)

Product Details

Sequence:

FPDGTKEMDI VFHEHQDNGT QDDSLTLASV NTDFAFSLYK KLALKNPDTN IVFSPLSISA
ALALVSLGAK GKTMEEILEG LKFNLTETPE ADIHQGFGNL LQSLSQPEDQ DQINIGNAMF
IEKDLQILAE FHEKTRALYQ TEAFTADFQQ PTEAKNLIND YVSNQTQGMI KELISELDER
TLMVLVNYIY FKGKWKISFD PQDTFESEFY LDEKRSVKVP MMKMKLLTTR HFRDEELSCS
VLELKYTGNA SALLILPDQG RMQQVEASLQ PETLRKWRKT LFPSQIEELN LPKFSIASNY
RLEEDVLPEM GIKEVFTEQA DLSGITETKK LSVSQVVHKA VLDVAETGTE AAAATGVIGG
IRKAILPAVH FNRPFLFVIY HTSAQSILFM AKVNNPK

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Serpina3k Protein (raised in E. Coli) purified by multi-step, protein-specific process to ensure crystallization grade.

• State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in bacterial culture:

- In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility:	0.22 µm filtered
Endotoxin Level:	Endotoxin has not been removed. Please contact us if you require endotoxin removal.
Grade:	Crystallography grade

Target Details

Target:	SERPINA3K
Alternative Name:	Serpina3k (SERPINA3K Products)

Target Details

Background:	Contrapsin inhibits trypsin-like proteases. {EC0:0000269 PubMed:6214866, EC0:0000269 PubMed:6224776}.
Molecular Weight:	45.7 kDa Including tag.
UniProt:	P07759

Application Details

Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies
	as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee
	though.
Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the
	recombinant protein with the default tag will be insoluble our protein lab may suggest a higher
	molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible
	options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)



Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process