

Datasheet for ABIN3133866

**Fibrillarin Protein (FBL) (AA 1-327) (His tag)**[Go to Product page](#)**1** Image

## Overview

|                               |  |
|-------------------------------|--|
| Quantity:                     | 1 mg   |
| Target:                       | Fibrillarin (FBL)  |
| Protein Characteristics:      | AA 1-327   |
| Origin:                       | Mouse  |
| Source:                       | Escherichia coli (E. coli)   |
| Protein Type:                 | Recombinant  |
| Purification tag / Conjugate: | This Fibrillarin protein is labelled with His tag.                   |
| Application:                  | SDS-PAGE (SDS), Western Blotting (WB), ELISA, Crystallization (Crys) |

## Product Details

|           |  |
|-----------|--|
| Sequence: | MKPGFSPRGG GFGGRRGGFGD RGGRRGGRRGG RGGFGGGRGG FGGGGRGRGG GGGGFRGRGG<br>GGGRGGGFQS GGNRRGGGR GKGKRNQSGK NVMVEPHRHE GVFCIRGKED ALVTKNLVPG<br>ESVYGEKRVS ISEGDDKIEY RAWNPFRSKL AAAILGGVDQ IHIKPGAKVL YLGAASGTTV<br>SHVSDIVGPD GLVYAVEFSH RSGRDLINLA KKRTNIIPVI EDARHPHKYR MLIAMVDVIF<br>ADVAQPDQTR IVALNAHTFL RGGGHFVISI KANCIDSTAS AEAVFASEVK KMQQENMKPQ<br>EQLTLEPYER DHAVVGVYR PPPKVKN |
|-----------|--|

**Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.**

|                  |   |
|------------------|---|
| Characteristics: | <ul style="list-style-type: none"><li>• Made in Germany - from design to production - by highly experienced protein experts.</li><li>• Mouse Fbl Protein (raised in E. Coli) purified by multi-step, protein-specific process to ensure crystallization grade.</li><li>• State-of-the-art algorithm used for plasmid design (Gene synthesis).</li></ul> |
|------------------|---|

## Product Details

---

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receipt of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the ExPASy's protparam tool to determine the absorption coefficient of each protein.

|                  |  |
|------------------|--|
| Purification:    | Two step purification of proteins expressed in bacterial culture:<br><ol style="list-style-type: none"><li>1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.</li><li>2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.</li></ol> |
| Purity:          | >95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.   |
| Sterility:       | 0.22 µm filtered   |
| Endotoxin Level: | Endotoxin has not been removed. Please contact us if you require endotoxin removal.  |
| Grade:           | Crystallography grade  |

## Target Details

---

|                   |  |
|-------------------|--|
| Target:           | Fibrillarin (FBL)  |
| Alternative Name: | Fbl ( <a href="#">FBL Products</a> )   |
| Background:       | S-adenosyl-L-methionine-dependent methyltransferase that has the ability to methylate both |

## Target Details

RNAs and proteins. Involved in pre-rRNA processing by catalyzing the site-specific 2'-hydroxyl methylation of ribose moieties in pre-ribosomal RNA. Site specificity is provided by a guide RNA that base pairs with the substrate. Methylation occurs at a characteristic distance from the sequence involved in base pairing with the guide RNA. Also acts as a protein methyltransferase by mediating methylation of 'Gln-105' of histone H2A (H2AQ104me), a modification that impairs binding of the FACT complex and is specifically present at 35S ribosomal DNA locus (By similarity). {ECO:0000250}.

Molecular Weight: 35.3 kDa Including tag.

UniProt: [P35550](#)

Pathways: [Ribonucleoside Biosynthetic Process](#)

## Application Details

Application Notes: In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.

Comment: Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.

Restrictions: For Research Use only

## Handling

Format: Liquid

Buffer: 100 mM NaCl, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.

Handling Advice: Avoid repeated freeze-thaw cycles.

Storage: -80 °C

Storage Comment: Store at -80°C.

Expiry Date: Unlimited (if stored properly)



**Image 1.** „Crystallography Grade“ protein due to multi-step, protein-specific purification process