

Datasheet for ABIN3134610

TMPRSS9 Protein (AA 53-1065) (His tag)[Go to Product page](#)**1** Image

Overview

Quantity:	1 mg
Target:	TMPRSS9
Protein Characteristics:	AA 53-1065
Origin:	Mouse
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This TMPRSS9 protein is labelled with His tag.
Application:	Crystallization (Crys), ELISA, SDS-PAGE (SDS), Western Blotting (WB)

Product Details

Sequence:	SAQGAHVEHI AELHGIRFTS SLQQENSDFY RLLTHALQTL LHFLLRALQP LSLKQEADIL QKGIQARLQE QGLSLAAYGT IVSVELTGRC EGPVTERDLK SGHCPGNVFS CQNGQCVSKE NPECDDRVC SDESDEAQCD CGWQPAWRSA GRIVGGVEAA PGEFPWQVSL RENHEHFCGA TIIGARWLVS AAHCFNEFQD PAQWAAQAGS VHLSGSEASA VRTRVLRIAK HPAYDADTAD FDVAVLELAR PLPFGRYVQP ACLPAATHVF PPGKKCLISG WGYLKEDFLV KPEVLQKATV ELLDQSLCSS LYGHSLTDRM VCAGYLDGKV DSCQGDSCGGP LVCEEPSGRF FLAGIVSWG GCAEARRPGV YTRVTRLRDW ILEVTSADM PVVPTATPAP ATPSTPWPTS PESWAPNTFA KPTAAPSPVP LHPSTTAKPQ ECGARPAMDK PTRIVGGISA VSGEVPWQAS LKEGPRHFCG ATVVGDRWLL SAAHCFNHTK VEQVQAH LGT VSL LGVGGSP VKLGLRRVAL HPRYNPGILD FDVALLEL AQ PLVFNKYIQ VCLPLAIHKF PVGRKCMISG WGNMQEGNAT KPDILQKASV GIIEQKMC GA LYNFSLTDRM LCAGFLEGRV DSCQGDSCGGP LACEETPGVF YLAGIVSWG GCAQAKKPGV YARITRLKDW ILKAMSSDPS SMARPHTSST RLIPSEPPKT TAAGLIPEA
-----------	--

TTSRLATPRA TIRVTTRPLN TTLSARSTTT RGQTAAPSAP GTTIHSHLPD CGLAPPGALT
RIVGGSAASL GEWPWQVSLW LRRREHRCGA VLVAERWLLS AAHCFDIYGD PMQWAAFLGT
PFLSSTEGQL ERVARIYRHP FYNIYTLDYD VALLELAGPV RRSRLVRPIC LPGPARPPDG
ARCVITGWGS LREGGSMARQ LQKAAVRVLS EQTCRRFYPV QISSRMLCAG FPQGGVDSCS
GDAGGPLACR EPSGQWVLTG VTSWGYGCGR PHFPGVYTRV AAVLGWIGQN IQE

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany - from design to production - by highly experienced protein experts.
- Mouse Tmprss9 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receipt of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the ExPASy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
2. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Product Details

Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility:	0.22 µm filtered
Endotoxin Level:	Protein is endotoxin free.
Grade:	Crystallography grade

Target Details

Target:	TMPRSS9
Alternative Name:	Tmprss9 (TMPRSS9 Products)
Background:	Serase-1 and serase-2 are serine proteases that hydrolyze the peptides N-t-Boc-Gln-Ala-Arg-AMC and N-t-Boc-Gln-Gly-Arg-AMC. In contrast, N-t-Boc-Ala-Phe-Lys-AMC and N-t-Boc-Ala-Pro-Ala-AMC are not significantly hydrolyzed (By similarity). {ECO:0000250}.
Molecular Weight:	110.2 kDa Including tag.
UniProt:	P69525

Application Details

Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only

Handling

Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.

Handling

Expiry Date: Unlimited (if stored properly)

Images



Image 1. „Crystallography Grade“ protein due to multi-step, protein-specific purification process