
# antibodies .- online.com





# MERTK Protein (AA 19-994) (rho-1D4 tag)



**Image** 



Go to Product page

#### Overview

Quantity:	1 mg
Target:	MERTK
Protein Characteristics:	AA 19-994
Origin:	Mouse
Source:	Insect Cells
Protein Type:	Recombinant
Purification tag / Conjugate:	This MERTK protein is labelled with rho-1D4 tag.
Application:	Western Blotting (WB), SDS-PAGE (SDS), ELISA, Crystallization (Crys)

#### **Product Details**

Sequence:

GGTAEKWEET ELDQLFSGPL PGRLPVNHRP FSAPHSSRDQ LPPPQTGRSH PAHTAAPQVT
STASKLLPPV AFNHTIGHIV LSEHKNVKFN CSINIPNTYQ ETAGISWWKD GKELLGAHHS
ITQFYPDEEG VSIIALFSIA SVQRSDNGSY FCKMKVNNRE IVSDPIYVEV QGLPYFIKQP
ESVNVTRNTA FNLTCQAVGP PEPVNIFWVQ NSSRVNEKPE RSPSVLTVPG LTETAVFSCE
AHNDKGLTVS KGVHINIKVI PSPPTEVHIL NSTAHSILVS WVPGFDGYSP LQNCSIQVKE
ADRLSNGSVM VFNTSASPHL YEIQQLQALA NYSIAVSCRN EIGWSAVSPW ILASTTEGAP
SVAPLNITVF LNESNNILDI RWTKPPIKRQ DGELVGYRIS HVWESAGTYK ELSEEVSQNG
SWAQIPVQIH NATCTVRIAA ITKGGIGPFS EPVNIIIPEH SKVDYAPSST PAPGNTDSMF
IILGCFCGFI LIGLILCISL ALRRRVQETK FGGAFSEEDS QLVVNYRAKK SFCRRAIELT LQSLGVSEEL
QNKLEDVVID RNLLVLGKVL GEGEFGSVME GNLKQEDGTS QKVAVKTMKL DNFSQREIEE
FLSEAACMKD FNHPNVIRLL GVCIELSSQG IPKPMVILPF MKYGDLHTFL LYSRLNTGPK
YIHLQTLLKF MMDIAQGMEY LSNRNFLHRD LAARNCMLRD DMTVCVADFG LSKKIYSGDY

YRQGRIAKMP VKWIAIESLA DRVYTSKSDV WAFGVTMWEI TTRGMTPYPG VQNHEMYDYL LHGHRLKQPE DCLDELYDIM YSCWSADPLD RPTFSVLRLQ LEKLSESLPD AQDKESIIYI NTQLLESCEG IANGPSLTGL DMNIDPDSII ASCTPGAAVS VVTAEVHENN LREERYILNG GNEEWEDVSS TPFAAVTPEK DGVLPEDRLT KNGVSWSHHS TLPLGSPSPD ELLFVDDSLE DSEVLM

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

#### Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Mertk Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

#### Purification:

Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:

- 1. Membrane proteins are fractioned by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.
- 2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.
- 3. Protein containing fractions of the best purification are subjected to second purification step


### **Product Details**

	through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.
Purity:	>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Sterility:	0.22 µm filtered
Endotoxin Level:	Protein is endotoxin-free.
Grade:	Crystallography grade
Target Details	
Target:	MERTK
Alternative Name:	Mertk (MERTK Products)
Background:	Receptor tyrosine kinase that transduces signals from the extracellular matrix into the
	cytoplasm by binding to several ligands including LGALS3, TUB, TULP1 or GAS6. Regulates
	many physiological processes including cell survival, migration, differentiation, and
	phagocytosis of apoptotic cells (efferocytosis). Ligand binding at the cell surface induces
	autophosphorylation of MERTK on its intracellular domain that provides docking sites for
	downstream signaling molecules. Following activation by ligand, interacts with GRB2 or PLCG2
	and induces phosphorylation of MAPK1, MAPK2, FAK/PTK2 or RAC1. MERTK signaling plays a
	role in various processes such as macrophage clearance of apoptotic cells, platelet
	aggregation, cytoskeleton reorganization and engulfment. Functions in the retinal pigment
	epithelium (RPE) as a regulator of rod outer segments fragments phagocytosis. Plays also an
	important role in inhibition of Toll-like receptors (TLRs)-mediated innate immune response by
	activating STAT1, which selectively induces production of suppressors of cytokine signaling
	SOCS1 and SOCS3. {ECO:0000269 PubMed:19117932, ECO:0000269 PubMed:20363878}.
Molecular Weight:	109.4 kDa Including tag.
UniProt:	Q60805
Pathways:	RTK Signaling
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies
	as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee

## **Application Details**

Comment:	Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)

## Images



**Image 1.** "Crystallography Grade" protein due to multi-step, protein-specific purification process