antibodies .- online.com







GLDN Protein (AA 39-549) (His tag)



Image



Overview

| Quantity: | 1 mg |
|-------------------------------|--|
| Target: | GLDN |
| Protein Characteristics: | AA 39-549 |
| Origin: | Mouse |
| Source: | Insect Cells |
| Protein Type: | Recombinant |
| Purification tag / Conjugate: | This GLDN protein is labelled with His tag. |
| Application: | ELISA, Western Blotting (WB), Crystallization (Crys), SDS-PAGE (SDS) |

Product Details

Sequence:

COWRGLSAAL RALEAORGRE OREDSALRAF LAELSRAPGR VPEPSODPMS AARNKRSHNG EPASHIRAES QDMMMMMTYS MVPIRVMIDL CNSTQGICLT GPPGPPGPPG AGGLPGHNGS DGQPGLQGPK GEKGAIGKRG KMGLPGATGN PGEKGEKGDA GELGLPGNEG PPGQKGDKGD KGDVSNDVLL TGAKGDQGPP GPPGPPGPPG PPGSRRSKGP RPPNVFNSQC PGETCVIPND DTLVGRADEK ANERHSPQTE SMITSIGNPA QVLKVRETFG TWMRESANKS DDRIWVTEHF SGIMVKEFKD LPALLNSSFT LLHLPHYFHG CGHAVYNNSL YYHKGGSNTI VRFEFGKETP QTLKLENALY FDRKYLFANS KTYFNIAVDE KGIWIIYASS VDGSSILVAQ LDERTFSVTQ HINTTYPKSK AGNAFIARGI LYVTDTKDTR VTFAFDLLGG KQINANFDFR MSQSVLAMLS YNMRDQHLYS WEDGHLMLYP VQFLSAASSQ R

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany from design to production by highly experienced protein experts.
- Mouse Gldn Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in baculovirus infected SF9 insect cells:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Protein is endotoxin free.

Grade:

Crystallography grade

Target Details

| Target: | GLDN |
|--|---|
| Alternative Name: | Gldn (GLDN Products) |
| Background: | Ligand for NRCAM and NFASC/neurofascin that plays a role in the formation and maintenance of the nodes of Ranvier on myelinated axons. Mediates interaction between Schwann cell microvilli and axons via its interactions with NRCAM and NFASC (PubMed:20188654). Nodes of Ranvier contain clustered sodium channels that are crucial for the saltatory propagation of action potentials along myelinated axons. During development, nodes of Ranvier are formed by the fusion of two heminodes. Required for normal clustering of sodium channels at heminodes, not required for the formation of mature nodes with normal sodium channel clusters (PubMed:20188654). Required, together with NRCAM, for maintaining NFASC and sodium channel clusters at mature nodes of Ranvier (PubMed:24719088). {ECO:0000250 UniProtKB:Q80WL1, ECO:0000269 PubMed:20188654, ECO:0000269 PubMed:24719088}. |
| Molecular Weight: | 56.3 kDa Including tag. |
| UniProt: | Q8BMF8 |
| | |
| Application Details | |
| Application Details Application Notes: | In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. |
| | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee |
| Application Notes: | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible |
| Application Notes: Comment: | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. |
| Application Notes: Comment: Restrictions: | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. |
| Application Notes: Comment: Restrictions: Handling | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. For Research Use only |
| Application Notes: Comment: Restrictions: Handling Format: | as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. For Research Use only Liquid |

Handling

| Storage Comment: | Store at -80°C. |
|------------------|--------------------------------|
| Expiry Date: | Unlimited (if stored properly) |

Images




Image 1. "Crystallography Grade" protein due to multi-step, protein-specific purification process