# antibodies .- online.com





# ACAP1 Protein (AA 1-740) (Strep Tag)



## Overview

Quantity:	1 mg
Target:	ACAP1
Protein Characteristics:	AA 1-740
Origin:	Mouse
Source:	Tobacco (Nicotiana tabacum)
Protein Type:	Recombinant
Purification tag / Conjugate:	This ACAP1 protein is labelled with Strep Tag.
Application:	ELISA, Western Blotting (WB), SDS-PAGE (SDS)

# **Product Details**

Sequence:

MTVKLDFEEC LKDSPRFRAS IELVETEVSE LETRLEKLLK LGSCLLESGQ HYLAAGRAFV
VGICDLARLG PPEPMMAECL EKFTVSLNHK LDSHAELLDA TQHTLQQQIQ TLVKEGLRGF
REARRDFWRG AENLEAALTH NAEVPRRRVQ EAEEAGTALR TARAGYRSRA LDYALQVNVI
EDKRKFDIME FVLRLVEAQA TYFQQGHEEL NRLAQYRKEL GTQLHNLVLN SARQKRDMEQ
RHVLLKQKEL GGEEPEPSLK EGPSGLVMEG HLFKRASNAF KTWSRRWFTI QNNQLVYQKK
YKDPVTVVVD DLRLCTVKLC PDSERRFCFE VVSTSKSCFL QADSERLLQL WVSAVQSSIA
SAFSQAHLEN SPRGPGQVSG YHAPGSAATL ACGGAARGRE SGGVGQVAAQ VQSVDGNAQC
CDCREPAPEW ASINLGVTLC IQCSGIHRSL GVHFSKVRSL TLDSWEPELV KLMCELGNVI
INQIYEARVE AMAVKKPGPS CSRQEKEAWI HAKYVEKKFL TKLPEIRGRR GGRGPPRGHP
PVPPKPPIRP HSGIVRSKSE CPSDDMGSLH PGALLFQAAG HPPSLPTMAD ALAHGADVNW
VNVGQGNATP LIRATAANSL LACEFLLQNG ANVNQADSAG RGPLHHATIL GHTGLACLFL
KRGADLGARD TEGRDPLTIA METTNADIVT LLRLAKMREA EAAQGQAGDE TYLDIFRDFS

#### LMASDDPEKL SRRSHDLHTL

Sequence without tag. The proposed Strep-Tag is based on experience s with the expression system, a different complexity of the protein could make another tag necessary. In case you have a special request, please contact us.

#### Characteristics:

# Key Benefits:

- · Made in Germany from design to production by highly experienced protein experts.
- Protein expressed with ALiCE® and purified by multi-step, protein-specific process to ensure correct folding and modification.
- These proteins are normally active (enzymatically functional) as our customers have reported (not tested by us and not guaranteed).
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a **made-to-order protein** and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our **made-to-order proteins** in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

# Expression System:

- ALiCE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v.. This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require posttranslational modifications.
- During lysate production, the cell wall and other cellular components that are not required for
  protein production are removed, leaving only the protein production machinery and the
  mitochondria to drive the reaction. During our lysate completion steps, the additional
  components needed for protein production (amino acids, cofactors, etc.) are added to
  produce something that functions like a cell, but without the constraints of a living system all that's needed is the DNA that codes for the desired protein!

## Concentration:

- · The concentration of our recombinant proteins is measured using the absorbance at 280nm.
- The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.
- We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in Almost Living Cell-Free Expression System

	(ALiCE®):
	<ol> <li>In a first purification step, the protein is purified from the cleared cell lysate using StrepTag capture material. Eluate fractions are analyzed by SDS-PAGE.</li> <li>Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.</li> </ol>
Purity:	≥ 80 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.
Endotoxin Level:	Low Endotoxin less than 1 EU/mg (< 0.1 ng/mg)
Target Details	
Target:	ACAP1
Alternative Name:	Acap1 (ACAP1 Products)
Background:	Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (Centaurin-beta-1) (Cnt-b1), FUNCTION: GTPase-activating protein (GAP) for ADP ribosylation factor 6 (ARF6) required for clathrin-dependent export of proteins from recycling endosomes to trans-Golgi network and cell surface. Required for regulated export of ITGB1 from recycling endosomes to the cell surface and ITGB1-dependent cell migration (By similarity). {ECO:0000250, ECO:0000269 PubMed:17664335}.
Molecular Weight:	81.7 kDa
UniProt:	Q8K2H4
Application Details	
Application Notes:	In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though.
Comment:	ALICE®, our Almost Living Cell-Free Expression System is based on a lysate obtained from Nicotiana tabacum c.v This contains all the protein expression machinery needed to produce even the most difficult-to-express proteins, including those that require post-translational modifications.  During lysate production, the cell wall and other cellular components that are not required for protein production are removed, leaving only the protein production machinery and the mitochondria to drive the reaction. During our lysate completion steps, the additional components needed for protein production (amino acids, cofactors, etc.) are added to produce

# **Application Details**

	something that functions like a cell, but without the constraints of a living system - all that's needed is the DNA that codes for the desired protein!
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	The buffer composition is at the discretion of the manufacturer. If you have a special request, please contact us.
Handling Advice:	Avoid repeated freeze-thaw cycles.
Storage:	-80 °C
Storage Comment:	Store at -80°C.
Expiry Date:	Unlimited (if stored properly)