

Datasheet for ABIN3137126

TMPRSS6 Protein (AA 1-811) (rho-1D4 tag)[Go to Product page](#)**1** Image

Overview

| | |
|-------------------------------|--|
| Quantity: | 1 mg |
| Target: | TMPRSS6 |
| Protein Characteristics: | AA 1-811 |
| Origin: | Mouse |
| Source: | Insect Cells |
| Protein Type: | Recombinant |
| Purification tag / Conjugate: | This TMPRSS6 protein is labelled with rho-1D4 tag. |
| Application: | Western Blotting (WB), SDS-PAGE (SDS), ELISA, Crystallization (Crys) |

Product Details

| | |
|-----------|--|
| Sequence: | <p>MPRCFQLPCS TRMPTTEVPQ AADGQGDAGD GEEAAEPEGK FKPPKNTKRK NRDYVRFTPL LLVLAALVSA GVMLWYFLGY KAEVTVSQVY SGSLRVLNRH FSQDLGRRES IAFRSESAKA QKMLQELVAS TRLGTYYNSS SVYSFGEGPL TCFFWFILDI PEYQRLTLSP EVVRELLVDE LLSNSSTLAS YKTEYEVDPE GLVILEASVN DIVVLNSTLG CYRYSYVNPQ QVLPLKGPDQ QTTSCWLHLQ GPEDLMIKVR LEWTRVDCRD RVAMYDAAGP LEKRLITSVY GCSRQEPVME VLASGSVMAS VWKKGMHSYY DPFLSVKSV AFQDCQVNLT LEGRDLDQGF LRTPYYPSYY SPSTHCSWHL TVPSLDYGLA LWFDAYALRR QKYNRLCTQG QWMIQNRRLC GFRTLQPYAE RIPMVASDGV TINFTSQISL TGPGVQVYYS LYNQSDPCPG EFLCSVNGLC VPACDGIKDC PNGLDERNVC CRAMFQCQED STCISLPRVC DRQPDCLNGS DEEQCQEGVP CGTFTFCQED RSCVKKPNPE CDGQSDCRDG SDEQHCDCGL QGLSSRIVGG TVSSEGEWPW QASLQIRGRH ICGGALIADR WVITAACHCFQ EDSMASPKLW TVFLGKMRQN SRWPGEVSFK VSRLFLHPYH EEDSHDYDVA LLQLDHPVYV SATVRPVCLP ARSHFFEPGQ HCWITGWGAQ REGGPVSNTL</p> |
|-----------|--|

QKVDVQLVPQ DLCSEAYRYQ VSPRMLCAGY RKGKKDACQG DSGGPLVCRE PSGRWFLAGL
VSWGLGCGRP NFFGVYTRVT RVINWIIQVL T

Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a special request, please contact us.

Characteristics:

- Made in Germany - from design to production - by highly experienced protein experts.
- Mouse Tmprss6 Protein (raised in Insect Cells) purified by multi-step, protein-specific process to ensure crystallization grade.
- State-of-the-art algorithm used for plasmid design (Gene synthesis).

This protein is a made to order protein and will be made for the first time for your order. Our experts in the lab will ensure that you receive a correctly folded protein.

The big advantage of ordering our made-to-order proteins in comparison to ordering custom made proteins from other companies is that there is no financial obligation in case the protein cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receipt of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered.

The concentration of our recombinant proteins is measured using the absorbance at 280nm.

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the ExPASy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Three step purification of membrane proteins expressed in baculovirus infected SF9 insect cells:

1. Membrane proteins are fractionated by ultracentrifugation and subsequently solubilized with different detergents (detergent screen). Samples are analyzed by Western blot.
2. The best performing detergent is used for solubilization and the proteins are purified via their rho1D4 tag via two rho1D4 antibody columns: one DTT resistant, the other one not. Eluate fractions are analyzed by Western blot.
3. Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatograph. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Product Details

| | |
|------------------|----------------------------|
| Sterility: | 0.22 µm filtered |
| Endotoxin Level: | Protein is endotoxin-free. |
| Grade: | Crystallography grade |

Target Details

| | |
|-------------------|--|
| Target: | TMPRSS6 |
| Alternative Name: | Tmprss6 (TMPRSS6 Products) |
| Background: | Serine protease which hydrolyzes a range of proteins including type I collagen, fibronectin and fibrinogen. Can also activate urokinase-type plasminogen activator with low efficiency. May play a specialized role in matrix remodeling processes in liver. Through the cleavage of HFE2, a regulator of the expression of the iron absorption-regulating hormone hepcidin/HAMP, plays a role in iron homeostasis. {ECO:0000269 PubMed:18451267}. |
| Molecular Weight: | 92.2 kDa Including tag. |
| UniProt: | Q9DBI0 |
| Pathways: | Transition Metal Ion Homeostasis |

Application Details

| | |
|--------------------|---|
| Application Notes: | In addition to the applications listed above we expect the protein to work for functional studies as well. As the protein has not been tested for functional studies yet we cannot offer a guarantee though. |
| Comment: | Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. |
| Restrictions: | For Research Use only |

Handling

| | |
|------------------|--|
| Format: | Liquid |
| Buffer: | 100 mM NaCl, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer. |
| Handling Advice: | Avoid repeated freeze-thaw cycles. |

Handling

Storage: -80 °C

Storage Comment: Store at -80°C.

Expiry Date: Unlimited (if stored properly)

Images



Image 1. „Crystallography Grade“ protein due to multi-step, protein-specific purification process