

Datasheet for ABIN3146684

SLC22A5 Protein (AA 42-142) (His tag)



Go to Product page

_					
	W	0	rv	10	W

Quantity:	2 mg		
Target:	SLC22A5		
Protein Characteristics:	AA 42-142		
Origin:	Mouse		
Source:	Escherichia coli (E. coli)		
Protein Type:	Recombinant		
Purification tag / Conjugate:	This SLC22A5 protein is labelled with His tag.		
Application:	ELISA, Western Blotting (WB), SDS-PAGE (SDS), Crystallization (Crys)		
Product Details			
Sequence:	LAGTPEHRCL VPHTVNLSSA WRNHSIPLET KDGRQVPQKC RRYRLATIAN FSELGLEPGR		
	DVDLEQLEQE SCLDGWEYDK DVFLSTIVTE WDLVCKDDWK A		
	Sequence without tag. Tag location is at the discretion of the manufacturer. If you have a		
	special request, please contact us.		
Characteristics:	 Made in Germany - from design to production - by highly experienced protein experts. Mouse Slc22a5 Protein (raised in E. Coli) purified by multi-step, protein-specific process to ensure crystallization grade. State-of-the-art algorithm used for plasmid design (Gene synthesis). 		
	This protein is a made to order protein and will be made for the first time for your order. Our		
	experts in the lab will ensure that you receive a correctly folded protein.		
	The big advantage of ordering our made-to-order proteins in comparison to ordering custom		
	ang annum go an annum go an annum pantum pantum go annum go annum go annum go annum go annum go annum go a		

cannot be expressed or purified.

In the unlikely event that the protein cannot be expressed or purified we do not charge anything (other companies might charge you for any performed steps in the expression process for custom-made proteins, e.g. fees might apply for the expression plasmid, the first expression experiments or purification optimization).

When you order this made-to-order protein you will only pay upon receival of the correctly folded protein. With no financial risk on your end you can rest assured that our experienced protein experts will do everything to make sure that you receive the protein you ordered. The concentration of our recombinant proteins is measured using the absorbance at 280nm. The protein's absorbance will be measured in several dilutions and is measured against its

The protein's absorbance will be measured in several dilutions and is measured against its specific reference buffer.

The concentration of the protein is calculated using its specific absorption coefficient. We use the Expasy's protparam tool to determine the absorption coefficient of each protein.

Purification:

Two step purification of proteins expressed in bacterial culture:

- 1. In a first purification step, the protein is purified from the cleared cell lysate using three different His-tag capture materials: high yield, EDTA resistant, or DTT resistant. Eluate fractions are analyzed by SDS-PAGE.
- Protein containing fractions of the best purification are subjected to second purification step through size exclusion chromatography. Eluate fractions are analyzed by SDS-PAGE and Western blot.

Purity:

>95 % as determined by SDS PAGE, Size Exclusion Chromatography and Western Blot.

Sterility:

0.22 µm filtered

Endotoxin Level:

Endotoxin has not been removed. Please contact us if you require endotoxin removal.

Grade:

Crystallography grade

Target Details

Target:	SLC22A5	
Alternative Name:	Slc22a5 (SLC22A5 Products)	
Background:	Sodium-ion dependent, high affinity carnitine transporter. Involved in the active cellular uptake of carnitine. Transports one sodium ion with one molecule of carnitine. Also transports organic cations such as tetraethylammonium (TEA) without the involvement of sodium. Also relative uptake activity ratio of carnitine to TEA is 11.3.	
Molecular Weight:	12.6 kDa Including tag.	

Target Details UniProt: Q9Z0E8 **Application Details** In addition to the applications listed above we expect the protein to work for functional studies Application Notes: as well. As the protein has not been tested for functional studies yet we cannot offer a gurantee though. Comment: Protein has not been tested for activity yet. In cases in which it is highly likely that the recombinant protein with the default tag will be insoluble our protein lab may suggest a higher molecular weight tag (e.g. GST-tag) instead to increase solubility. We will discuss all possible options with you in detail to assure that you receive your protein of interest. Restrictions: For Research Use only Handling

Format:	Liquid	
Buffer:	100 mM NaCL, 20 mM Hepes, 10% glycerol. pH value is at the discretion of the manufacturer.	
Handling Advice:	Avoid repeated freeze-thaw cycles.	
Storage:	-80 °C	
Storage Comment:	Store at -80°C.	
Expiry Date:	Unlimited (if stored properly)	