# ANTIBODIES ONLINE

# Datasheet for ABIN366593 **Pyridinoline ELISA Kit**

1 Image

1 Publication



#### Overview

Quantity:	96 tests
Target:	Pyridinoline (PYD)
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	12.5-800 ng/mL
Minimum Detection Limit:	12.5 ng/mL
Application:	ELISA

## Product Details

Purpose:	For the quantitative determination of endogenic human pyridinoline (PYD) concentrations in serum, plasma, urine.
Sample Type:	Serum, Plasma, Urine
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of human PYD.
Cross-Reactivity (Details):	Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between the target antigen and all analogues for other species. Therefore, cross reaction may still exist.
Sensitivity:	3.12 ng/mL
Components:	Assay plate (12 × 8 coated Microwells)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/5 | Product datasheet for ABIN366593 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

- Standard (freeze dried)
- Biotin-antibody (100 × concentrate)
- HRP-avidin (100 × concentrate)
- Biotin-antibody Diluent
- HRP-avidin Diluent
- Sample Diluent
- Wash Buffer (25 × concentrate)
- TMB Substrate
- Stop Solution
- Adhesive Strip (for 96 wells)
- Instruction manual

### Target Details


Target:	Pyridinoline (PYD)
Alternative Name:	Pyridinoline crosslinks (PYD) (PYD Products)
Target Type:	Amino Acid
Application Details	
Application Notes:	<ul> <li>The supplier is only responsible for the kit itself, but not for the samples consumed during the assay. The user should calculate the possible amount of the samples used in the whole test. Please reserve sufficient samples in advance.</li> <li>Samples to be used within 5 days may be stored at 2-8°C, otherwise samples must be stored at -20°C (≤ 1 month) or -80°C (≤ 2 months) to avoid loss of bioactivity and contamination.</li> <li>Grossly hemolyzed samples are not suitable for use in this assay.</li> <li>If the samples are not indicated in the manual, a preliminary experiment to determine the validity of the kit is necessary.</li> <li>Please predict the concentration before assaying. If values for these are not within the range of the standard curve, users must determine the optimal sample dilutions for their particular experiments.</li> <li>Tissue or cell extraction samples prepared by chemical lysis buffer may cause unexpected ELISA results due to the impacts of certain chemicals.</li> <li>Owing to the possibility of mismatching between antigens from another resource and antibodies used in this supplier's kits (e.g., antibody targets conformational epitope rather than linear epitope), some native or recombinant proteins from other manufacturers may not be recognized by this supplier's products.</li> <li>Influenced by factors including cell viability, cell number and cell sampling time, samples from cell culture supernatant may not be recognized by the kit.</li> </ul>
	• Fresh samples without long time storage are recommended for the test. Otherwise, protein degradation and denaturalization may occur in those samples and finally lead to wrong

# Application Details

	results.
Comment:	Detection wavelength: 450 nm
	Information on standard material:
	Depending on the antigen to be detected, standards can be either native or recombinant
	protein. The recombinant proteins are being expressed in CHO cells in most cases. Please
	inquire for more information. The formulation of auxiliary material in the standard is considered
	proprietary information, however it does not contain any poisonous substance. Proclin 300
	(1:3000) is used as preservative.
	Information on reagents:
	In most cases the stop solution provided is 1 N H2SO4. The formulation of wash solution is
	proprietary information. None of the components contain (sodium) azide, thimerosal, 2-
	mercaptoethanol (2-ME) or any other poisonous materials. For the sandwich method kits, the
	sample diluent, antibody diluent, enzyme diluent and standard all contain BSA.
	Information on antibodies:
	The antibodies provided in different kits vary in regards to clonality and host. Some antibodies
	are affinity purified, some are Protein A
Sample Volume:	100 µL
Assay Time:	1 - 4.5 h
Plate:	Pre-coated
Protocol:	This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody
	specific for PYD has been pre-coated onto a microplate. Standards and samples are pipetted
	into the wells and any PYD present is bound by the immobilized antibody. After removing any
	unbound substances, a biotin-conjugated antibody specific for PYD is added to the wells. After
	washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a
	wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells
	and color develops in proportion to the amount of PYD bound in the initial step. The color
	development is stopped and the intensity of the color is measured.
Assay Precision:	Intra-assay precision (precision within an assay): Three samples of known concentration were
	tested twenty times on one plate to assess precision.
	Inter-assay precision (precision between assays): Three samples of known concentration were

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/5 | Product datasheet for ABIN366593 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Handling         recaution of Use:       The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing protection when using this material.         landling Advice:       • The kit should not be used beyond the expiration date on the kit label.         • Do not mix or substitute reagents with those from other lots or sources.       • If samples generate values higher than the highest standard, dilute the samples with Sampl Diluent and repeat the assay.         • Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.         • This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the limmunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone	Application Details	
<ul> <li>Inter-assay: CV% less than 10%</li> <li>estrictions: For Research Use only</li> <li>Handling</li> <li>recaution of Use: The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing protection when using this material.</li> <li>landling Advice: The kit should not be used beyond the expiration date on the kit label.</li> <li>Do not mix or substitute reagents with those from other lots or sources.</li> <li>If samples generate values higher than the highest standard, dilute the samples with Sampl Diluent and repeat the assay.</li> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.</li> <li>This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the immunoassay, the possibility of interference cannot be excluded.</li> <li>torage: 4 °C/-20 °C</li> <li>torage Comment: For unopened kit. All the reagents should be kept according to the labels on vials.</li> <li>xpiry Date: 6 months</li> <li>cublications</li> <li>roduct cited in: Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: According to the labels on vials.</li> </ul>		tested in twenty assays to assess precision.
estrictions:       For Research Use only         Handling       The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing protection when using this material.         Iandling Advice:       The Kit should not be used beyond the expiration date on the kit label.         Iandling Advice:       The kit should not be used beyond the expiration date on the kit label.         Iandling Advice:       The kit should not be used beyond the expiration date on the kit label.         Iandling Advice:       The kit should not be used beyond the expiration date on the kit label.         Iandling Advice:       The kit should not be used beyond the expiration date on the kit label.         Iandling Advice:       The kit should not be used beyond the expiration in the sor sources.         If samples generate values higher than the highest standard, dilute the samples with Sample Diluent and repeat the assay.       Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation time/temperature and kit age can cause variation in binding.         This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The eff		Intra-assay: CV% less than 8%
Handling         recaution of Use:       The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing protection when using this material.         landling Advice:       • The kit should not be used beyond the expiration date on the kit label.         andling Advice:       • The kit should not be used beyond the expiration date on the kit label.         bo not mix or substitute reagents with those from other lots or sources.       • If samples generate values higher than the highest standard, dilute the samples with Sample Diluent and repeat the assay.         Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.       • This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xplry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ageneration in postmenopausal osteoporotic women: a randomized clinical trial." in: Ageneration in postmenopausal osteoporotic women: a randomized clinical trial." in: Ageneration in postmenopausal osteoporotic women: a randomized clinical trial." in: Ageneration in postmenopausal osteoporotic w		Inter-assay: CV% less than 10%
recaution of Use:       The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing protection when using this material.         landling Advice:       • The kit should not be used beyond the expiration date on the kit label.         • Do not mix or substitute reagents with those from other lots or sources.       • If samples generate values higher than the highest standard, dilute the samples with Sample Diluent and repeat the assay.         • Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.         • This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomised clinical	Restrictions:	For Research Use only
protection when using this material.         landling Advice:         • The kit should not be used beyond the expiration date on the kit label.         • Do not mix or substitute reagents with those from other lots or sources.         • If samples generate values higher than the highest standard, dilute the samples with Sampl Diluent and repeat the assay.         • Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.         • This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postmenopausal osteoporotic women: a randomized clinical trial." in: Actional postremotion postremotichapausi postmenopausal postmenopausal postmeno	Handling	
<ul> <li>The kit should not be used beyond the expiration date on the kit label.</li> <li>Do not mix or substitute reagents with those from other lots or sources.</li> <li>If samples generate values higher than the highest standard, dilute the samples with Sampl Diluent and repeat the assay.</li> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.</li> <li>This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.</li> <li>torage Comment:</li> <li>For unopened kit: All the reagents should be kept according to the labels on vials.</li> <li>a months</li> </ul>	Precaution of Use:	The Stop Solution provided with this kit is an acid solution. Wear eye, hand, face and clothing
<ul> <li>Do not mix or substitute reagents with those from other lots or sources.</li> <li>If samples generate values higher than the highest standard, dilute the samples with Sampl Diluent and repeat the assay.</li> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation time/temperature and kit age can cause variation in binding.</li> <li>This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.</li> <li>torage: 4 °C/-20 °C</li> <li>torage Comment: For unopened kit: All the reagents should be kept according to the labels on vials.</li> <li>xpiry Date: 6 months</li> <li>coulciations</li> </ul>		protection when using this material.
<ul> <li>If samples generate values higher than the highest standard, dilute the samples with Sample Diluent and repeat the assay.</li> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubatio time/temperature and kit age can cause variation in binding.</li> <li>This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.</li> <li>torage: 4 °C/-20 °C</li> <li>torage Comment: For unopened kit: All the reagents should be kept according to the labels on vials.</li> <li>publications</li> </ul>	Handling Advice:	The kit should not be used beyond the expiration date on the kit label.
Diluent and repeat the assay.Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation time/temperature and kit age can cause variation in binding.This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.torage:4 °C/-20 °Ctorage Comment:For unopened kit: All the reagents should be kept according to the labels on vials.xpiry Date:6 monthsPublicationsSalari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Age		<ul> <li>Do not mix or substitute reagents with those from other lots or sources.</li> </ul>
time/temperature and kit age can cause variation in binding.         This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Age		
<ul> <li>This assay is designed to eliminate interference by soluble receptors, binding proteins and other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.</li> <li>torage: 4 °C/-20 °C</li> <li>torage Comment: For unopened kit: All the reagents should be kept according to the labels on vials.</li> <li>xpiry Date: 6 months</li> <li>Publications</li> <li>roduct cited in: Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag</li> </ul>		• Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation
other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.torage:4 °C/-20 °Ctorage Comment:For unopened kit: All the reagents should be kept according to the labels on vials.xpiry Date:6 monthsPublicationsSalari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag		time/temperature and kit age can cause variation in binding.
Immunoassay, the possibility of interference cannot be excluded.         torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag		This assay is designed to eliminate interference by soluble receptors, binding proteins and
torage:       4 °C/-20 °C         torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag		other factors present in biological samples. Until all factors have been tested in the
torage Comment:       For unopened kit: All the reagents should be kept according to the labels on vials.         xpiry Date:       6 months         Publications       5         roduct cited in:       Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag		Immunoassay, the possibility of interference cannot be excluded.
xpiry Date: 6 months Publications roduct cited in: Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag	Storage:	4 °C/-20 °C
Publications roduct cited in: Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag	Storage Comment:	For unopened kit: All the reagents should be kept according to the labels on vials.
roduct cited in: Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag	Expiry Date:	6 months
biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Ag	Publications	
	Product cited in:	Salari Sharif, Asalforoush, Ameri, Larijani, Abdollahi: "The effect of n-3 fatty acids on bone
(Dordrecht, Netherlands), Vol. 32, Issue 2, pp. 179-86, (2010) (PubMed).		biomarkers in Iranian postmenopausal osteoporotic women: a randomized clinical trial." in: Age
		(Dordrecht, Netherlands), Vol. 32, Issue 2, pp. 179-86, (2010) (PubMed).



#### ELISA

Image 1. Typical standard curve

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 5/5 | Product datasheet for ABIN366593 | 07/26/2024 | Copyright antibodies-online. All rights reserved.