

# Datasheet for ABIN367377 **PTX3 ELISA Kit**

1 Image

7 Publications



#### Overview

Quantity:	96 tests
Target:	PTX3
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156-10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

## Product Details

Purpose:	For the quantitative determination of human pentraxin 3 (PTX3) concentrations in serum, cell culture supernates, tissue homogenates.
Sample Type:	Serum, Cell Culture Supernatant, Tissue Homogenate, Plasma
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of human PTX3.
Cross-Reactivity (Details):	Limited by current skills and knowledge, it is impossible for us to complete the cross-reactivity detection between the target antigen and all analogues for other species. Therefore, cross reaction may still exist.
Sensitivity:	0.039 ng/mL
Components:	Assay plate (12 × 8 coated Microwells)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/6 | Product datasheet for ABIN367377 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

- Standard (freeze dried)
- Biotin-antibody (100 × concentrate)
- HRP-avidin (100 × concentrate)
- Biotin-antibody Diluent
- HRP-avidin Diluent
- Sample Diluent
- Wash Buffer (25 × concentrate)
- TMB Substrate
- Stop Solution
- Adhesive Strip (for 96 wells)
- Instruction manual

## Target Details

Target:	PTX3
Abstract:	PTX3 Products
Background:	Synonyms: TNFAIP5, TSG-14, pentaxin-related gene, rapidly induced by IL-1 beta, tumor necrosis factor, alpha-induced protein 5 pentaxin-related protein PTX3 pentraxin 3 pentraxin- 3 tumor necrosis factor, alph
HGNC:	9006
UniProt:	P26022

## Application Details

Application Notes:	• The supplier is only responsible for the kit itself, but not for the samples consumed during the
	assay. The user should calculate the possible amount of the samples used in the whole test.
	Please reserve sufficient samples in advance.
	• Samples to be used within 5 days may be stored at 2-8°C, otherwise samples must be stored
	at -20°C ( $\leq$ 1 month) or -80°C ( $\leq$ 2 months) to avoid loss of bioactivity and contamination.
	<ul> <li>Grossly hemolyzed samples are not suitable for use in this assay.</li> </ul>
	If the samples are not indicated in the manual, a preliminary experiment to determine the
	validity of the kit is necessary.
	• Please predict the concentration before assaying. If values for these are not within the range
	of the standard curve, users must determine the optimal sample dilutions for their particular
	experiments.
	Tissue or cell extraction samples prepared by chemical lysis buffer may cause unexpected
	ELISA results due to the impacts of certain chemicals.
	Owing to the possibility of mismatching between antigens from another resource and
	antibodies used in this supplier's kits (e.g., antibody targets conformational epitope rather
	than linear epitope), some native or recombinant proteins from other manufacturers may not

be recognized by this supplier's products.

- Influenced by factors including cell viability, cell number and cell sampling time, samples from cell culture supernatant may not be recognized by the kit.
- Fresh samples without long time storage are recommended for the test. Otherwise, protein degradation and denaturalization may occur in those samples and finally lead to wrong results.

Comment:

Detection wavelength: 450 nm

Information on standard material:

Depending on the antigen to be detected, standards can be either native or recombinant protein. The recombinant proteins are being expressed in CHO cells in most cases. Please inquire for more information. The formulation of auxiliary material in the standard is considered proprietary information, however it does not contain any poisonous substance. Proclin 300 (1:3000) is used as preservative.

Information on reagents:

In most cases the stop solution provided is 1 N H2SO4. The formulation of wash solution is proprietary information. None of the components contain (sodium) azide, thimerosal, 2-mercaptoethanol (2-ME) or any other poisonous materials. For the sandwich method kits, the sample diluent, antibody diluent, enzyme diluent and standard all contain BSA.

Information on antibodies:

The antibodies provided in different kits vary in regards to clonality and host. Some antibodies are affinity purified, some are Protein A

Sample Volume:	100 µL
Assay Time:	1 - 4.5 h
Plate:	Pre-coated
Protocol:	This assay employs the quantitative sandwich enzyme immunoassay technique. Antibody
	specific for PTX3 has been pre-coated onto a microplate. Standards and samples are pipetted
	into the wells and any PTX3 present is bound by the immobilized antibody. After removing any
	unbound substances, a biotin-conjugated antibody specific for PTX3 is added to the wells. After
	washing, avidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a
	wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells
	and color develops in proportion to the amount of PTX3 bound in the initial step. The color
	development is stopped and the intensity of the color is measured.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/6 | Product datasheet for ABIN367377 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

Restrictions:	For Research Use only
	<ul> <li>Intra-assay: CV% less than 8%</li> <li>Inter-assay: CV% less than 10%</li> </ul>
	tested in twenty assays to assess precision.
	Inter-assay precision (precision between assays): Three samples of known concentration were
Assay Precision:	Intra-assay precision (precision within an assay): Three samples of known concentration were tested twenty times on one plate to assess precision.
	contaminated water or container for reagent preparation will influence detection result.
	<ul> <li>It is recommended to use distilled water to prepare reagents and samples. Using</li> </ul>
	pipetting, use small volumes and ensure that pipettors are calibrated. It is recommended to suck more than 10µL when pipetting.
	gently until the crystals have completely dissolved. To minimize imprecision caused by
	Please carefully reconstitute Standards according to the instruction. Avoid foaming and mix
	Making serial dilution in the wells directly is not permitted.
	Prepare fresh standard for each assay. Use within 4 hours and discard after use.
	<ul> <li>Bring all reagents to room temperature (18-25°C) before use for 30 min.</li> </ul>
	directly in the Diluent vials provided in the kit.
	<ul> <li>Kindly use graduated containers to prepare the reagent. Please don't prepare the reagent</li> </ul>
	Note:
	serves as the high standard. Sample Diluent serves as the zero standard (Ong/mL).
	dilution series. Mix each tube thoroughly before the next transfer. The undiluted Standard
	making dilutions. Pipette 250µL of Sample Diluent into each tube. Use the stock solution to produce a 2-fold
	and allow the standard to sit for a minimum of 15 minutes with gentle agitation prior to
	reconstitution produces a stock solution. Mix the standard to ensure complete reconstitution
	Reconstitute the Standard with 1ml of Sample Diluent. Do not substitute other diluents. This
	• Standard - Centrifuge the standard vial at 6000-10000rpm for 30s.
	Concentrate (25×) into deionized or distilled water to prepare 500mL of Wash Buffer (1×).
	and mix gently until the crystals have completely dissolved. Dilute 20mL of Wash Buffer
	<ul> <li>of HRP-avidin Diluent.</li> <li>Wash Buffer (1×) - If crystals have formed in the concentrate, warm up to room temperature</li> </ul>
	HRP-avidin requires a 100-fold dilution. The suggested dilution is 10 $\mu$ L of HRP-avidin + 990 $\mu$ l
Reagent Preparation:	HRP-avidin (1×) - Centrifuge the vial before opening.
	+ 990µL of Biotin-antibody Diluent.
	<ul> <li>Biotin-antibody (1×) - Centrifuge the vial before opening.</li> <li>Biotin-antibody requires a 100-fold dilution. The suggested dilution is 10µL of Biotin-antibody</li> </ul>

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 4/6 | Product datasheet for ABIN367377 | 07/26/2024 | Copyright antibodies-online. All rights reserved.


	protection when using this material.
Handling Advice:	<ul> <li>The kit should not be used beyond the expiration date on the kit label.</li> </ul>
	<ul> <li>Do not mix or substitute reagents with those from other lots or sources.</li> <li>If samples generate values higher than the highest standard, dilute the samples with Sample</li> </ul>
	Diluent and repeat the assay.
	<ul> <li>Any variation in Sample Diluent, operator, pipetting technique, washing technique, incubation</li> </ul>
	time/temperature and kit age can cause variation in binding.
	This assay is designed to eliminate interference by soluble receptors, binding proteins and
	other factors present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.
	ininditoassay, the possibility of interference cannot be excluded.
Storage:	4 °C/-20 °C
Storage Comment:	For unopened kit: All the reagents should be kept according to the labels on vials.
Expiry Date:	6 months
Publications	
Product cited in:	Li, Liu, Wang, Shen, Du, Liu, Chen, Hu: "High apolipoprotein M serum levels correlate with
	chronic obstructive pulmonary disease." in: Lipids in health and disease, Vol. 15, pp. 59, (2016
	(PubMed).
	Rauf, Soejono, Partadiredja: "Effects of treadmill exercise training on cerebellar estrogen and
	estrogen receptors, serum estrogen, and motor coordination performance of ovariectomized
	rats." in: Iranian journal of basic medical sciences, Vol. 18, Issue 6, pp. 587-92, (2015) (PubMe
	).
	Pérez, Palma, Burgos, Jeldres, Espinoza, Peñailillo: "The acute phase response induced by
	Escherichia coli lipopolysaccharide modifies the pharmacokinetics and metabolism of
	florfenicol in rabbits." in: Journal of veterinary pharmacology and therapeutics, (2015) (
	PubMed).
	Chen, Chen, Liu, Sun, Mao: "Effect of Chromium Supplementation on Element Distribution in a
	Mouse Model of Polycystic Ovary Syndrome." in: Biological trace element research, (2015) (
	PubMed).
	Haybar, Assareh, Ghotbi, Torabizadeh, Bozorgmanesh: "Incremental diagnostic value of
	circulating pentraxin in patients with intermediate risk of coronary artery disease." in: Heart

Page 5/6 | Product datasheet for ABIN367377 | 07/26/2024 | Copyright antibodies-online. All rights reserved.

(British Cardiac Society), Vol. 99, Issue 9, pp. 640-8, (2013) (PubMed).

There are more publications referencing this product on: Product page

#### Images



### ELISA

Image 1. Typical standard curve

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 6/6 | Product datasheet for ABIN367377 | 07/26/2024 | Copyright antibodies-online. All rights reserved.