

Datasheet for ABIN377038
Propidium Iodide (PI)



[Go to Product page](#)

Overview

Quantity: 100 tests

Application: Flow Cytometry (FACS)

Product Details

Characteristics: ApoScreen® Propidium Iodide (PI)

Purification: Purified

Application Details

Application Notes:

- **Applications:** FC - Quality tested
- **Working Dilutions:** Flow Cytometry Propidium iodide 10 L/10⁶ cells For flow cytometry, the suggested use of these reagents is in a final volume of 100 L

Restrictions: For Research Use only

Handling

Buffer: 100 tests in 1.0 mL of ApoScreen® Annexin V Binding Buffer and a stabilizing agent

Handling Advice: Each reagent is stable for the period shown on the bottle label if stored as directed.

Storage: 4 °C

Storage Comment: Store at 2-8°C