



[Go to Product page](#)

Datasheet for ABIN377039
Annexin V Binding Buffer 10X

Overview

Quantity:	50 mL
Application:	Flow Cytometry (FACS)

Product Details

Characteristics:	ApoScreen® Annexin V Binding Buffer 10X
------------------	---

Application Details

Application Notes:	<ul style="list-style-type: none">• Applications: FC• Working Dilutions: Flow Cytometry Dilute 1 part 10X buffer with 9 parts dH2O just prior to use
--------------------	---

Restrictions:	For Research Use only
---------------	-----------------------

Handling

Format:	Liquid
Buffer:	ApoScreen® Annexin V Binding Buffer 10X is supplied as 50 mL.
Handling Advice:	Wash cells twice in cold PBS, remove the PBS from the cell pellet after the second wash. Each reagent is stable for the period shown on the bottle label if stored as directed.
Storage:	4 °C
Storage Comment:	Store at 2-8°C