

# Datasheet for ABIN411262

## **BMP2 ELISA Kit**





### Overview

Quantity:	96 tests
Target:	BMP2
Binding Specificity:	AA 283-396
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	62.5-4000 pg/mL
Minimum Detection Limit:	62.5 pg/mL
Application:	ELISA

### **Product Details**

Purpose:	Sandwich High Sensitivity ELISA kit for Quantitative Detection of Mouse BMP-2
Brand:	PicoKine™
Sample Type:	Osseous Tissue, Cell Culture Supernatant, Serum
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Immunogen:	Expression system for standard: CHO Immunogen sequence: Q283-R396
Specificity:	Expression system for standard: CHO Immunogen sequence: Q283-R396
Cross-Reactivity (Details):	There is no detectable cross-reactivity with other relevant proteins.

Product Details	
Predicted Reactivity:	Bovine,Canine,Horse,Monkey,Rabbit
Sensitivity:	<2pg/mL
Material not included:	Microplate reader in standard size. Automated plate washer. Adjustable pipettes and pipette tips. Multichannel pipettes are recommended in the condition of large amount of samples in the detection. Clean tubes and Eppendorf tubes. Washing buffer (neutral PBS or TBS). Preparation of 0.01M TBS: Add 1.2g Tris, 8.5g Nacl
Target Details	
Target:	BMP2
Alternative Name:	BMP2 (BMP2 Products)
Background:	Protein Function: Induces cartilage and bone formation.  Background: Bone morphogenetic protein-2(BMP-2, BMP 2A) belongs to the transforming growth factor beta(TGF-beta) superfamily. It is thought to be involved in cartilage and bone formation during embryogenesis, but may have additional functions in morphogenesis as implied by its expression in various organs and embryonic tissues of mice. BMP-2 has been identified as a candidate mediator of retinoid activity. BMP-2 protein induces medulloblastoma cell apoptosis, whereas the BMP-2 antagonist noggin blocks both retinoid and BMP-2-induced apoptosis. BMP-2 also induces p38 mitogen-activated protein kinase(MAPK), which is necessary for BMP-2 and retinoid-induced apoptosis. Bone morphogenic proteins(BMPs) are known to promote osteogenesis, and clinical trials are currently underway to evaluate the ability of certain BMPs to promote fracture-healing and spinal fusion. The standard product used in this kit is recombinant mouse BMP-2, constituting dimer by two chains of 114 amino acids with the molecular mass of 26KDa.  Synonyms: Bone morphogenetic protein 2,BMP-2,Bone morphogenetic protein 2A,BMP-

2A,Bmp2,Bmp-2,

Full Gene Name: Bone morphogenetic protein 2

Cellular Localisation: Secreted.

Gene ID: 12156

UniProt: P21274

Pathways: Regulation of Hormone Metabolic Process, Regulation of Hormone Biosynthetic Process,

Regulation of Muscle Cell Differentiation, Growth Factor Binding, Positive Regulation of fat Cell

Differentiation

## **Application Details**

Application Notes:	Before using Kit, spin tubes and bring down all components to bottom of tube. Duplicate well
	assay was recommended for both standard and sample testing.
Plate:	Pre-coated
Protocol:	mouse BMP-2 ELISA Kit was based on standard sandwich enzyme-linked immune-sorbent
	assay technology. A monoclonal antibody from rat specific for BMP-2 has been precoated ont
	96-well plates. Standards(CHO, Q283-R396) and test samples are added to the wells, a
	biotinylated detection polyclonal antibody from goat specific for BMP-2 is added subsequently
	and then followed by washing with PBS or TBS buffer. Avidin-Biotin-Peroxidase Complex was
	added and unbound conjugates were washed away with PBS or TBS buffer. HRP substrate
	TMB was used to visualize HRP enzymatic reaction. TMB was catalyzed by HRP to produce a
	blue color product that changed into yellow after adding acidic stop solution. The density of
	yellow is proportional to the mouse BMP-2 amount of sample captured in plate.
Assay Procedure:	Aliquot 0.1 mL per well of the 4000pg/mL, 2000pg/mL,1000pg/mL, 500pg/mL, 250pg/mL,
	125pg/mL, 62.5pg/mL mouse BMP-2 standard solutions into the precoated 96-well plate. Add
	0.1 mL of the sample diluent buffer into the control well (Zero well). Add 0.1 mL of each
	properly diluted sample of mouse bone tissue, cell culture supernates or serum to each empty
	well. See "Sample Dilution Guideline" above for details. It is recommended that each mouse
	BMP-2 standard solution and each sample be measured in duplicate.
Assay Precision:	• Sample 1: n=16, Mean(pg/ml): 314, Standard deviation: 11.3, CV(%): 3.6
	<ul> <li>Sample 2: n=16, Mean(pg/ml): 1002, Standard deviation: 39.1, CV(%): 3.9</li> </ul>
	• Sample 3: n=16, Mean(pg/ml): 2263, Standard deviation: 88.3, CV(%): 3.9,
	• Sample 1: n=24, Mean(pg/ml): 334, Standard deviation: 24, CV(%): 7.2
	<ul> <li>Sample 2: n=24, Mean(pg/ml): 1059, Standard deviation: 74.1, CV(%): 7</li> <li>Sample 3: n=24, Mean(pg/ml): 2360, Standard deviation: 177, CV(%): 7.5</li> </ul>
	Sumple 6.11 24, Wean(pg/111). 2000, Standard deviation. 177, 6 v (%). 7.5
Restrictions:	For Research Use only
Handling	
Handling Advice:	Avoid multiple freeze-thaw cycles.
Storage:	-20 °C,4 °C
Storage Comment:	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles
Expiry Date:	12 months

# 

#### **ELISA**

Image 1. Mouse BMP-2 PicoKine ELISA Kit standard curve