

Datasheet for ABIN4972481

Reelin ELISA Kit



Overview

| Quantity: | 96 tests |
|--------------------------|---------------------------|
| Target: | Reelin (RELN) |
| Reactivity: | Mouse |
| Method Type: | Sandwich ELISA |
| Detection Range: | 15.625 pg/mL - 1000 pg/mL |
| Minimum Detection Limit: | 15.625 pg/mL |
| Application: | ELISA |

Product Details

| Purpose: | Mouse Reelin ELISA Kit is an ELISA Kit against Reelin. |
|--------------------|--------------------------------------------------------|
| Sample Type: | Plasma, Serum |
| Analytical Method: | Quantitative |
| Detection Method: | Colorimetric |
| Sensitivity: | 9.38 pg/mL |

Target Details

| Target: | Reelin (RELN) |
|-------------------|------------------------|
| Alternative Name: | Reelin (RELN Products) |
| Pathways: | Synaptic Membrane |

Application Details

| Application Notes: | Stability: The stability of the kit is determined by the rate of activity loss. The loss rate is less |
|--------------------|-------------------------------------------------------------------------------------------------------|
| | than 5 % within the expiration date under appropriate storage conditions. To minimize |
| | performance fluctuations, operation procedures and lab conditions should be strictly controlled. |
| | It is also strongly suggested that the whole assay is performed by the same user throughout. |
| | Recommended dilutions: Optimal dilutions/concentrations should be determined by the end |
| | user. |
| | Standard Form: Lyophilized |
| Comment: | The stability of the kit is determined by the rate of activity loss. The loss rate is less than 5% |
| | within the expiration date under appropriate storage conditions. To minimize performance |
| | fluctuations, operation procedures and lab conditions should be strictly controlled. It is also |
| | strongly suggested that the whole assay is performed by the same user throughout. |
| Plate: | Pre-coated |
| Restrictions: | For Research Use only |
| Handling | |
| Storage: | 4 °C/-20 °C |
| Storage Comment: | Upon receipt, store the kit according to the storage instruction in the kit's manual. |
| Expiry Date: | 6 months |