# antibodies -online.com





## anti-COX3 antibody (AA 25-130) (Alexa Fluor 680)



Go to Product page

( )	11/0	K\ /	iew	1
	$\cup$	'I V/I	$\square \vee \vee$	ı

Quantity:	100 μL	
Target:	COX3 (COX-3)	
Binding Specificity:	AA 25-130	
Reactivity:	Mouse, Rat	
Host:	Rabbit	
Clonality:	Polyclonal	
Conjugate:	This COX3 antibody is conjugated to Alexa Fluor 680	
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))	

### **Product Details**

Immunogen:	KLH conjugated synthetic peptide derived from human COX3
Isotype:	IgG
Cross-Reactivity:	Mouse, Rat
Predicted Reactivity:	Human,Dog,Cow,Pig,Horse,Chicken,Rabbit
Purification:	Purified by Protein A.

#### **Target Details**

Target:	COX3 (COX-3)
Alternative Name:	COX3 (COX-3 Products)

### **Target Details**

rarget Details		
Background:	Synonyms: Cytochrome c oxidase polypeptide III, Cytochrome c oxidase subunit 3, COIII, COX3, COX-3, COXIII, MT CO3, MTCO3, COX3_HUMAN.	
	Background: Cytochrome C Oxidase subunits I, II and III form the functional core of the enzyme	
	complex.Defects in Cytochrome C Oxidase subunits III are a cause of Leber hereditary optic	
	neuropathy (LHON) and cytochrome c oxidase deficiency (COX deficiency). Defects in	
	Cytochrome C Oxidase subunits III are also found in mitochondrial encephalomyopathy with	
	lactic acidosis and stroke-like episodes (MELAS) syndrome and recurrent myoglobinuria.	
Gene ID:	4514	
Application Details		
Application Notes:	IF(IHC-P) 1:50-200	
	IF(IHC-F) 1:50-200	
	IF(ICC) 1:50-200	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and	
	50 % Glycerol.	
Preservative:	ProClin	
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be	
	handled by trained staff only.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	
Expiry Date:	12 months	