antibodies -online.com







anti-PMP22 antibody (AA 101-160) (AbBy Fluor® 750)



Go to Product page

()	1/0	r\ /1	014	
()	ve	I V I	-v	V

Quantity:	100 μL
Target:	PMP22
Binding Specificity:	AA 101-160
Reactivity:	Rat, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This PMP22 antibody is conjugated to AbBy Fluor® 750
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human PMP-22
Isotype:	IgG
Cross-Reactivity:	Mouse, Rat
Predicted Reactivity:	Human
Purification:	Purified by Protein A.

Target Details

Target:	PMP22
Alternative Name:	PMP22 (PMP22 Products)

Target Details

Background:	Synonyms: DSS, HNPP, CMT1A, CMT1E, GAS-3, Sp110, HMSNIA, Peripheral myelin protein 22,	
	PMP-22, Growth arrest-specific protein 3, PMP22, GAS3	
	Background: PMP22 is a 22 kDa glycoprotein expressed in the compact myelin of the peripheral	
	nervous system. In the peripheral nervous system, PMP 22 is produced by myelinating	
	Schwann cells and is coexpressed with the genes for myelin basic protein (MBP) during nerve	
	development and regeneration. Alterations in the level of this protein cause several genetic	
	human diseases. If the protein is duplicated, patients develop Charcot Marie Tooth disease. If	
	one copy of the gene is deleted, they suffer from the inherited tendency to pressure palsies.	
Gene ID:	5376	
UniProt:	Q01453	
Application Details		
Application Notes:	IF(IHC-P) 1:50-200	
	IF(IHC-F) 1:50-200	
	IF(ICC) 1:50-200	
Restrictions:	For Research Use only	
Handling		
Format:	Liquid	
Concentration:	1 μg/μL	
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and	
	50 % Glycerol.	
Preservative:	ProClin	
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be	
	handled by trained staff only.	
Storage:	-20 °C	
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.	
Expiry Date:	12 months	