# antibodies -online.com





## anti-Supervillin antibody (AA 51-150) (Alexa Fluor 750)



Overview	
Quantity:	100 μL
Target:	Supervillin (SVIL)
Binding Specificity:	AA 51-150
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Supervillin antibody is conjugated to Alexa Fluor 750
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

#### **Product Details**

Immunogen:	KLH conjugated synthetic peptide derived from human Supervillin	
Isotype:	IgG	
Cross-Reactivity:	Mouse	
Predicted Reactivity:	Human,Rat,Dog,Cow,Sheep,Pig,Horse,Chicken	
Purification:	Purified by Protein A.	

#### **Target Details**

Target:	Supervillin (SVIL)	
Alternative Name:	Supervillin (SVIL Products)	

### **Target Details**

Background:	Synonyms: Supervillin, Archvillin, p205/p250, SVIL		
	Background: Forms a high-affinity link between the actin cytoskeleton and the membrane.		
	Isoform 1 (archvillin) is among the first costameric proteins to assemble during myogenesis		
	and it contributes to myogenic membrane structure and differentiation. Appears to be involved in myosin II assembly. May modulate myosin II regulation through MLCK during cell spreading,		
	involved in modulation of focal adhesions. Supervillin-mediated down-regulation of focal		
	adhesions involves binding to TRIP6. Plays a role in cytokinesis through KIF14 interaction (By		
	similarity).		
Gene ID:	6840		
UniProt:	095425		
Application Details			
Application Notes:	IF(IHC-P) 1:50-200		
	IF(IHC-F) 1:50-200		
	IF(ICC) 1:50-200		
Restrictions:	For Research Use only		
Handling			
Format:	Liquid		
Concentration:	1 μg/μL		
Buffer:	Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and		
	50 % Glycerol.		
Preservative:	ProClin		
Precaution of Use:	This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be		
	handled by trained staff only.		
Storage:	-20 °C		
Storage Comment:	Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.		
Expiry Date:	12 months		