

Datasheet for ABIN5012212

**anti-HAUS7 antibody (AA 41-140) (AbBy Fluor® 750)**[Go to Product page](#)

## Overview

Quantity:	100 µL
Target:	HAUS7
Binding Specificity:	AA 41-140
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This HAUS7 antibody is conjugated to AbBy Fluor® 750
Application:	Western Blotting (WB), Immunofluorescence (Cultured Cells) (IF (cc)), Immunofluorescence (Paraffin-embedded Sections) (IF (p))

## Product Details

Immunogen:	KLH conjugated synthetic peptide derived from human UCHL5IP
Isotype:	IgG
Predicted Reactivity:	Human, Mouse, Rat, Dog, Cow, Sheep, Horse, Rabbit
Purification:	Purified by Protein A.

## Target Details

Target:	HAUS7
Alternative Name:	UCHL5IP ( <a href="#">HAUS7 Products</a> )
Background:	Synonyms: 26S proteasome-associated UCH interacting protein 1, 26S proteasome-associated

## Target Details

---

UCH37-interacting protein 1, Expressed-Xq28STS protein, HAUS augmin-like complex subunit 7, Haus7, HAUS7\_HUMAN, UCH37 interacting protein 1, UCHL5 interacting protein, UCHL5-interacting protein, UCHL5IP, UIP1, X-linked protein STS1769.

Background: Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex. Tissue specificity: Detected in spleen, thymus, testis, ovary, small intestine and colon, with highest levels of expression in testis and ovary.

---

Gene ID: 55559

## Application Details

---

Application Notes: IF(IHC-P) 1:50-200  
IF(IHC-F) 1:50-200  
IF(ICC) 1:50-200

---

Restrictions: For Research Use only

## Handling

---

Format: Liquid

---

Concentration: 1 µg/µL

---

Buffer: Aqueous buffered solution containing 0.01M TBS (pH 7.4) with 1 % BSA, 0.03 % Proclin300 and 50 % Glycerol.

---

Preservative: ProClin

---

Precaution of Use: This product contains ProClin: a POISONOUS AND HAZARDOUS SUBSTANCE, which should be handled by trained staff only.

---

Storage: -20 °C

---

Storage Comment: Store at -20°C. Aliquot into multiple vials to avoid repeated freeze-thaw cycles.

---

Expiry Date: 12 months