# antibodies - online.com







# anti-Manic Fringe antibody (AA 56-309)





( )	ve	K\ /		A .
	$\cup$	1 V/	Щ.	V۷

Quantity:	100 μL
Target:	Manic Fringe (MFNG)
Binding Specificity:	AA 56-309
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This Manic Fringe antibody is un-conjugated
Application:	Western Blotting (WB), Immunohistochemistry (IHC), Immunocytochemistry (ICC)

# **Product Details**

Immunogen:	MFNG (Ile56-Tyr309)
Isotype:	IgG
Specificity:	The antibody is a rabbit polyclonal antibody raised against MFNG. It has been selected for its ability to recognize MFNG in immunohistochemical staining and western blotting.
Purification:	Antigen-specific affinity chromatography

# **Target Details**

Target:	Manic Fringe (MFNG)	
Alternative Name:	Beta-1,3-N-Acetylglucosaminyltransferase Manic Fringe (MFNG) (MFNG Products)	
Background: Alternative Names: MFNG O-Fucosylpeptide 3-Beta-N-Acetylglucosaminyltransferase		

Pathways:

Notch Signaling

# **Application Details**

### Application Notes:

Western blotting: 1:50-400 Immunocytochemistry in formalin fixed cells: 1:50-500
 Immunohistochemistry in formalin fixed frozen section: 1:50-500 Immunohistochemistry in paraffin section: 1:10-100 Enzyme-linked Immunosorbent Assay: 1:100-1:5000 Optimal working dilutions must be determined by end user.

### Comment:

The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37&degC for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.


Restrictions:

For Research Use only

# Handling


Format:	Liquid
Concentration:	Lot specific
Buffer:	PBS, pH 7.4, containing 0.02 % Sodium azide, 50 % glycerol.

# **Images**




## **Western Blotting**

**Image 1.** Figure. Western Blot; Sample: Recombinant protein.



# **Immunohistochemistry**

**Image 2.** DAB staining on IHC-P; Samples: Rat Spleen Tissue



# **Immunohistochemistry**

**Image 3.** Used in DAB staining on fromalin fixed paraffinembedded Kidney tissue