



[Go to Product page](#)

Datasheet for ABIN5514402  
**anti-HDA11 antibody (Middle Region)**

### Overview

Quantity:	100 µL
Target:	HDA11 (HDA-11)
Binding Specificity:	Middle Region
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Conjugate:	This HDA11 antibody is un-conjugated
Application:	Western Blotting (WB)

### Product Details

Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of Human HDA11
Sequence:	RGGGFCAYAD ITLAIKFLFE RVEGISRATI IDLDAHQNGNG HERDFMDDKR
Characteristics:	This is a rabbit polyclonal antibody against HDA11. It was validated on Western Blot.
Purification:	Affinity purified

### Target Details

Target:	HDA11 (HDA-11)
Alternative Name:	HDA11 ( <a href="#">HDA-11 Products</a> )
Background:	This gene encodes a class IV histone deacetylase. The encoded protein is localized to the nucleus and may be involved in regulating the expression of interleukin 10. Alternative splicing

## Target Details

---

results in multiple transcript variants.

Alias Symbols: HDAC11,

Protein Interaction Partner: ATP2A3, RBM12B, ZNF326, SLC25A42, RCOR2, WDR36, PAXBP1, CIRH1A, CLPB, NUP85, HDAC11, RNF219, C2orf47, HHIP, MCCC2, XAB2, WRNIP1, PPAN, YLPM1, UTP6, NUP133, RBM22, DDX27, TTC27, HJURP, WDR33, PPP6R3, RIF1, RBM28, ACIN1, RUNX1, CALU, CAD, BDH1, IKZF1, DNAJA

Protein Size: 347

---

Gene ID: 79885

---

NCBI Accession: [NP\\_079103](#)

---

UniProt: [Q96DB2](#)

## Application Details

---

Application Notes: Optimal working dilution should be determined by the investigator.

Restrictions: For Research Use only

## Handling

---

Format: Liquid

Buffer: Liquid. Purified antibody supplied in 1x PBS buffer with 0.09 % (w/v) sodium azide and 2 % sucrose.

Preservative: Sodium azide

Precaution of Use: This product contains Sodium azide: a POISONOUS AND HAZARDOUS SUBSTANCE which should be handled by trained staff only.

Storage: -20 °C

Storage Comment: For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.