

Datasheet for ABIN5564559
Arylsulfatase A ELISA Kit



[Go to Product page](#)

1 Image

Overview

Quantity: 96 tests

Target: Arylsulfatase A (ARSA)

Reactivity: Human

Method Type: Sandwich ELISA

Application: ELISA

Product Details

Brand: AssayMax™

Sample Type: Cell Culture Cells, Plasma, Serum

Analytical Method: Quantitative

Detection Method: Colorimetric

Material not included: Microplate reader capable of measuring absorbance at 405 nm. Pipettes (1-20 µL, 20-200 µL, and multiple channel). Deionized or distilled reagent grade water Incubator (37 °C)

Target Details

Target: Arylsulfatase A (ARSA)

Alternative Name: Arylsulfatase A (ASA) ([ARSA Products](#))

Gene ID: 410

UniProt: [P15289](#)

Application Details

Plate: Pre-coated

Calculation of Results:

- Calculate the mean value of the duplicate or triplicate readings for each standard and sample.
- To generate a standard curve, plot the graph using the standard concentrations on the x-axis and the corresponding mean 450 nm absorbance (OD) on the y-axis. The best-fit line can be determined by regression analysis using log-log or four-parameter logistic curve-fit.
- Determine the unknown sample concentration from the standard curve and multiply the value by the dilution factor.

Restrictions: For Research Use only

Handling

Handling Advice: This product is for Research Use Only and is Not For Use In Diagnostic Procedures. Prepare all reagents (working diluent buffer, wash buffer, standard, biotinylated antibody, and SP Conjugate) as instructed, prior to running the assay. Prepare all samples prior to running the assay. The dilution factors for the samples are suggested in this insert. However, the user should determine the optimal dilution factor. Spin down the SP Conjugate vial and the biotinylated antibody vial before opening and using contents. The Stop Solution is an acidic solution. The kit should not be used beyond the expiration date. 2

Storage: 4 °C, -20 °C

Images

Human Arylsulfatase A Standard Curve




Image 1.