

Datasheet for ABIN5651269

OAS2 ELISA Kit



Overview

Quantity:	96 tests
Target:	OAS2
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.781 ng/mL - 50 ng/mL
Minimum Detection Limit:	0.781 ng/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of 2',5'-Oligoadenylate Synthetase 2 (OAS2). No significant cross-reactivity or interference between 2',5'-Oligoadenylate Synthetase 2 (OAS2) and analogues was observed.
Sensitivity:	0.36 ng/mL

Target Details

Target:	OAS2
Alternative Name:	2',5'-Oligoadenylate Synthetase 2 (OAS2 Products)

Target Details

Restrictions:

Background:	Gene Name: 2',5'-Oligoadenylate Synthetase 2
	Gene Aliases: p69 OAS / p71 OAS
Gene ID:	4939
UniProt:	P29728
Pathways:	Hepatitis C
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to 2',5'-Oligoadenylate
	Synthetase 2 (OAS2). Standards or samples are then added to the appropriate microtiter plate
	wells with a biotin-conjugated antibody specific to 2',5'-Oligoadenylate Synthetase 2 (OAS2).
	Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and
	incubated. After TMB substrate solution is added, only those wells that contain 2',5'-
	Oligoadenylate Synthetase 2 (OAS2), biotin-conjugated antibody and enzyme-conjugated Avidir
	will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of
	sulphuric acid solution and the color change is measured spectrophotometrically at a
	wavelength of 450nm \pm 10nm. The concentration of 2',5'-Oligoadenylate Synthetase 2 (OAS2) in
	the samples is then determined by comparing the O.D. of the samples to the standard curve.
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	2',5'-Oligoadenylate Synthetase 2 (OAS2) were tested 20 times on one plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	2',5'-Oligoadenylate Synthetase 2 (OAS2) were tested on 3 different plates, 8 replicates in each
	plate. CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%
D	5 B 111 1

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn | International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com | Page 2/3 | Product datasheet for ABIN5651269 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

For Research Use only

Handling

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months