

Datasheet for ABIN5651877 Apelin 17 ELISA Kit



Overview

Quantity:	96 tests
Target:	Apelin 17 (AP17)
Reactivity:	Human
Method Type:	Competition ELISA
Detection Range:	24.69 pg/mL - 2000 pg/mL
Minimum Detection Limit:	24.69 pg/mL
Application:	ELISA

Product Details

Sample Type:	Cell Culture Supernatant, Cell Lysate, Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Apelin 17 (AP17). No significant cross-reactivity or interference between Apelin 17 (AP17) and analogues was observed.
Sensitivity:	9.14 pg/mL
Target Details	
Target:	Apelin 17 (AP17)
Alternative Name:	Apelin 17 (AP17 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5651877 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details

Background:

Gene Name: Apelin 17 Gene Aliases: APLN17

Application Details

Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	2 h
Plate:	Pre-coated
Protocol:	This assay employs the competitive inhibition enzyme immunoassay technique. A monoclonal
	antibody specific to Apelin 17 (AP17) has been pre-coated onto a microplate. A competitive
	inhibition reaction is launched between biotin labeled Apelin 17 (AP17) and unlabeled Apelin 17
	(AP17) (Standards or samples) with the pre-coated antibody specific to Apelin 17 (AP17). After
	incubation the unbound conjugate is washed off. Next, avidin conjugated to Horseradish
	Peroxidase (HRP) is added to each microplate well and incubated. The amount of bound HRP
	conjugate is reverse proportional to the concentration of Apelin 17 (AP17) in the sample. After
	addition of the substrate solution, the intensity of color developed is reverse proportional to the
	concentration of Apelin 17 (AP17) in the sample.
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	Apelin 17 (AP17) were tested 20 times on one plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Apelin 17 (AP17) were tested on 3 different plates, 8 replicates in each plate. $CV(\%)$ =
	SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5651877 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Handling

Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 3/3 | Product datasheet for ABIN5651877 | 07/25/2024 | Copyright antibodies-online. All rights reserved.