

# Datasheet for ABIN5652177

# **CD57 ELISA Kit**



Go to Product page

| ()     | 11/ | $\sim$ | r\ | ń | 0 | ۱۸/ | , |
|--------|-----|--------|----|---|---|-----|---|
| $\cup$ | V   | C      | rv | ı | ヒ | ٧V  |   |

| Quantity:                | 96 tests               |
|--------------------------|------------------------|
| Target:                  | CD57 (B3GAT1)          |
| Reactivity:              | Human                  |
| Method Type:             | Sandwich ELISA         |
| Detection Range:         | 0.156 ng/mL - 10 ng/mL |
| Minimum Detection Limit: | 0.156 ng/mL            |
| Application:             | ELISA                  |

## **Product Details**

| Sample Type:       | Cell Lysate, Tissue Homogenate   |
|--------------------|--|
| Analytical Method: | Quantitative   |
| Detection Method:  | Colorimetric   |
| Specificity:       | This assay has high sensitivity and excellent specificity for detection of Beta-1,3-Glucuronylt ransferase 1 (b3GAT1). No significant cross-reactivity or interference between Beta-1,3-Glucuronyltransferase 1 (b3GAT1) and analogues was observed. |
| Sensitivity:       | 0.053 ng/mL  |

# **Target Details**

| Target:           | CD57 (B3GAT1)                                      |
|-------------------|--|
| Alternative Name: | Beta-1,3-Glucuronyltransferase 1 (B3GAT1 Products) |

# Target Details

| Background:         | Gene Name: Beta-1,3-Glucuronyltransferase 1   |
|---------------------|---|
|                     | Gene Aliases: CD57, B3-GAT1, GLCATP, GlcAT-P, GlcUAT-P, HNK-1, LEU7, NK-1,                              |
|                     | Glucuronosyltransferase P, Galactosylgalactosylxylosylprotein 3-Beta-Glucuronosyltransferase            |
|                     | 1   |
| Pathways:           | Glycosaminoglycan Metabolic Process   |
| Application Details |   |
| Comment:            | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than |
|                     | 5 % within the expiration date under appropriate storage condition. To minimize extra influence         |
|                     | on the performance, operation procedures and lab conditions, especially room temperature, air           |
|                     | humidity, incubator temperature should be strictly controlled. It is also strongly suggested that       |
|                     | the whole assay is performed by the same operator from the beginning to the end.                        |
| Assay Time:         | 3 h   |
| Plate:              | Pre-coated  |
| Protocol:           | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate             |
|                     | provided in this kit has been pre-coated with an antibody specific to Beta-1,3-                         |
|                     | Glucuronyltransferase 1 (b3GAT1). Standards or samples are then added to the appropriate                |
|                     | microtiter plate wells with a biotin-conjugated antibody specific to Beta-1,3-                          |
|                     | Glucuronyltransferase 1 (b3GAT1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is            |
|                     | added to each microplate well and incubated. After TMB substrate solution is added, only those          |
|                     | wells that contain Beta-1,3-Glucuronyltransferase 1 (b3GAT1), biotin-conjugated antibody and            |
|                     | enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is               |
|                     | terminated by the addition of sulphuric acid solution and the color change is measured                  |
|                     | spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Beta-1,3-                  |
|                     | Glucuronyltransferase 1 (b3GAT1) in the samples is then determined by comparing the O.D. of             |
|                     | the samples to the standard curve.  |
| Assay Precision:    | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level            |
|                     | Beta-1,3-Glucuronyltransferase 1 (b3GAT1) were tested 20 times on one plate, respectively               |
|                     | Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level             |
|                     | Beta-1,3-Glucuronyltransferase 1 (b3GAT1) were tested on 3 different plates, 8 replicates in            |
|                     | each plate. CV(%) = SD/meanX100   |
|                     | Intra-Assay: CV<10%   |
|                     | Inter-Assay: CV<12%   |

# **Application Details**

| Restrictions:    | For Research Use only   |
|------------------|---|
| Handling         |   |
| Handling Advice: | The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and  |
|                  | specimen samples should be assayed in duplicate. Once the procedure has been started, all     |
|                  | steps should be completed without interruption.   |
| Storage:         | 4 °C,-20 °C   |
| Storage Comment: | -20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at |
|                  | 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant   |
|                  | pack. Minimize freeze/thaw cycles.  |
| Expiry Date:     | 4-8 months  |