

## Datasheet for ABIN5652240

## **BST1 ELISA Kit**



| Overview                 |   |
|--------------------------|---|
| Quantity:                | 96 tests  |
| Target:                  | BST1  |
| Reactivity:              | Human   |
| Method Type:             | Sandwich ELISA  |
| Detection Range:         | 78.12 pg/mL - 5000 pg/mL  |
| Minimum Detection Limit: | 78.12 pg/mL   |
| Application:             | ELISA   |
| Product Details          |   |
| Sample Type:             | Plasma, Serum   |
| Analytical Method:       | Quantitative  |
| Detection Method:        | Colorimetric  |
| Specificity:             | This assay has high sensitivity and excellent specificity for detection of Bone Marrow Stromal Cell Antigen 1 (BST1). No significant cross-reactivity or interference between Bone Marrow Stromal Cell Antigen 1 (BST1) and analogues was observed. |
| Sensitivity:             | 31 pg/mL  |
| Target Details           |   |
| Target:                  | BST1  |
| Alternative Name:        | Bone Marrow Stromal Cell Antigen 1 (BST1 Products)  |

## **Target Details**

| Background:         | Gene Name: Bone Marrow Stromal Cell Antigen 1  |  |  |
|---------------------|--|--|--|
|                     | Gene Aliases: CD157, ADP-Ribosyl Cyclase 2, ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase  |  |  |
|                     | 2  |  |  |
| Gene ID:            | 683  |  |  |
| UniProt:            | Q10588   |  |  |
| Application Details |  |  |  |
| Comment:            | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.   |  |  |
| Assay Time:         | 3 h  |  |  |
| Plate:              | Pre-coated   |  |  |
| Protocol:           | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Bone Marrow Stromal Antigen 1 (BST1). Standards or samples are then added to the appropriate microtiter plate with a biotin-conjugated antibody specific to Bone Marrow Stromal Cell Antigen 1 (BST1). Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Bone Marrow Stromal Cell Antigen 1 (BST1), biotin-conjugated antibody and enzyme-conjugated Avidin exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Bone Marrow Stromal Cell Antigen 1 (BST1) in the samples is then determined by comparing the 0.D. of the samples to the stacurve. |  |  |
| Assay Precision:    | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Bone Marrow Stromal Cell Antigen 1 (BST1) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Bone Marrow Stromal Cell Antigen 1 (BST1) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%  |  |  |

## **Application Details**

| Restrictions:    | For Research Use only   |  |  |
|------------------|---|--|--|
| Handling         |   |  |  |
| Handling Advice: | The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and  |  |  |
|                  | specimen samples should be assayed in duplicate. Once the procedure has been started, all     |  |  |
|                  | steps should be completed without interruption.   |  |  |
| Storage:         | 4 °C,-20 °C   |  |  |
| Storage Comment: | -20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at |  |  |
|                  | 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant   |  |  |
|                  | pack. Minimize freeze/thaw cycles.  |  |  |
| Expiry Date:     | 4-8 months  |  |  |