

## Datasheet for ABIN5652263

### **BMP3 ELISA Kit**



#### Overview

Quantity:	96 tests
Target:	BMP3
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	31.25 pg/mL - 2000 pg/mL
Minimum Detection Limit:	31.25 pg/mL
Application:	ELISA

#### **Product Details**

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Bone Morphogenetic Protein 3 (BMP3). No significant cross-reactivity or interference between Bone Morphogenetic Protein 3 (BMP3) and analogues was observed.
Sensitivity:	12.2 pg/mL

### **Target Details**

Target:	BMP3
Alternative Name:	Bone Morphogenetic Protein 3 (BMP3 Products)

# **Target Details** Gene Name: Bone Morphogenetic Protein 3 Background: Gene Aliases: BMP3A, BMP-3A, Osteogenic Gene ID: 651 UniProt: P12645 **Application Details** The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than Comment: 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. Assay Time: 3 h Plate: Pre-coated Protocol: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Bone Morphogenetic Protein 3 (BMP3). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Bone Morphogenetic Protein 3 (BMP3). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and

sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Bone Morphogenetic Protein 3 (BMP3) in the samples is then determined by comparing the O.D. of the samples to the standard curve.

Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Bone Morphogenetic Protein 3 (BMP3) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Bone Morphogenetic Protein 3 (BMP3) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10% Inter-Assay: CV<12%

incubated. After TMB substrate solution is added, only those wells that contain Bone

Morphogenetic Protein 3 (BMP3), biotin-conjugated antibody and enzyme-conjugated Avidin

will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn | International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com | Page 2/3 | Product datasheet for ABIN5652263 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

For Research Use only

Assay Precision:

Restrictions:

## Handling

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months