

Datasheet for ABIN5652342

BCAR1 ELISA Kit



Overview

Quantity:	96 tests
Target:	BCAR1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	78.12 pg/mL - 5000 pg/mL
Minimum Detection Limit:	78.12 pg/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Breast Cancer Anti- Estrogen Resistance 1 (BCAR1). No significant cross-reactivity or interference between Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) and analogues was observed.
Sensitivity:	27 pg/mL

Target Details

Target:	BCAR1
Alternative Name:	Breast Cancer Anti-Estrogen Resistance 1 (BCAR1 Products)

Target Details

Background:	Gene Name: Breast Cancer Anti-Estrogen Resistance 1	
	Gene Aliases: CAS, CRKAS, P130Cas, CASS1, Crk-Associated Substrate, Cas Scaffolding Protein	
	Family Member 1	
Gene ID:	9564	
UniProt:	P56945	
Pathways:	EGFR Signaling Pathway, Neurotrophin Signaling Pathway, CXCR4-mediated Signaling Events,	
	Platelet-derived growth Factor Receptor Signaling	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Breast Cancer Anti-	
	Estrogen Resistance 1 (BCAR1). Standards or samples are then added to the appropriate	
	microtiter plate wells with a biotin-conjugated antibody specific to Breast Cancer Anti-Estrogen	
	Resistance 1 (BCAR1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to	
	each microplate well and incubated. After TMB substrate solution is added, only those wells	
	that contain Breast Cancer Anti-Estrogen Resistance 1 (BCAR1), biotin-conjugated antibody and	
	enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is	
	terminated by the addition of sulphuric acid solution and the color change is measured	
	spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Breast Cancer	
	Anti-Estrogen Resistance 1 (BCAR1) in the samples is then determined by comparing the O.D.	
	of the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) were tested 20 times on one plate,	
	respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Breast Cancer Anti-Estrogen Resistance 1 (BCAR1) were tested on 3 different plates, 8	

Application Details

	replicates in each plate. CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at
	4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant
	pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months