

## Datasheet for ABIN5652353 BRD4 ELISA Kit



Overview

Quantity:	96 tests
Target:	BRD4
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

## Product Details

Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Bromodomain Containing Protein 4 (BRD4). No significant cross-reactivity or interference between Bromodomain Containing Protein 4 (BRD4) and analogues was observed.
Sensitivity:	0.059 ng/mL
Target Details	

Target:	BRD4	
Alternative Name:	Bromodomain Containing Protein 4 (BRD4 Products)	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5652353 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
Background:	Gene Name: Bromodomain Containing Protein 4	
	Gene Aliases: CAP, HUNKI, MCAP, HUNK1, Chromosome-Associated Protein	
Pathways:	Chromatin Binding, SARS-CoV-2 Protein Interactome	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Bromodomain Containing	
	Protein 4 (BRD4). Standards or samples are then added to the appropriate microtiter plate wells	
	with a biotin-conjugated antibody specific to Bromodomain Containing Protein 4 (BRD4). Next,	
	Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and	
	incubated. After TMB substrate solution is added, only those wells that contain Bromodomain	
	Containing Protein 4 (BRD4), biotin-conjugated antibody and enzyme-conjugated Avidin will	
	exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of	
	sulphuric acid solution and the color change is measured spectrophotometrically at a	
	wavelength of 450nm $\pm$ 10nm. The concentration of Bromodomain Containing Protein 4 (BRD4)	
	in the samples is then determined by comparing the O.D. of the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Bromodomain Containing Protein 4 (BRD4) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Bromodomain Containing Protein 4 (BRD4) were tested on 3 different plates, 8 replicates in	
	each plate. CV(%) = SD/meanX100	
	Intra-Assay: CV<10%	
	Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5652353 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months