

## Datasheet for ABIN5652496 CAMK2A ELISA Kit



Overview

Quantity:	96 tests
Target:	CAMK2A
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

## Product Details

Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2a). No significant cross-reactivity or interference between Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2a) and analogues was observed.
Sensitivity:	0.118 ng/mL

## Target Details

Target:	CAMK2A
Alternative Name:	Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2A Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5652496 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Background:	Gene Name: Calcium/Calmodulin Dependent Protein Kinase II Alpha
	Gene Aliases: CAMKA
Pathways:	WNT Signaling, Interferon-gamma Pathway, Myometrial Relaxation and Contraction
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, ai
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to Calcium/Calmodulin
	Dependent Protein Kinase II Alpha (CAMK2a). Standards or samples are then added to the
	appropriate microtiter plate wells with a biotin-conjugated antibody specific to
	Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2a). Next, Avidin conjugated to
	Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB
	substrate solution is added, only those wells that contain Calcium/Calmodulin Dependent
	Protein Kinase II Alpha (CAMK2a), biotin-conjugated antibody and enzyme-conjugated Avidin
	will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of
	sulphuric acid solution and the color change is measured spectrophotometrically at a
	wavelength of 450nm $\pm$ 10nm. The concentration of Calcium/Calmodulin Dependent Protein
	Kinase II Alpha (CAMK2a) in the samples is then determined by comparing the O.D. of the
	samples to the standard curve.
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2a) were tested 20 times on on
	plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Calcium/Calmodulin Dependent Protein Kinase II Alpha (CAMK2a) were tested on 3 different
	plates, 8 replicates in each plate. CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5652496 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Application Details	
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months