

Datasheet for ABIN5652937 CGA ELISA Kit



Overview

Quantity:	96 tests
Target:	CGA
Reactivity:	Human
Method Type:	Competition ELISA
Detection Range:	123.5 pg/mL - 10000 pg/mL
Minimum Detection Limit:	123.5 pg/mL
Application:	ELISA

Product Details

Sample Type:	Cell Culture Supernatant, Cell Lysate, Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Chorionic Gonadotropir Alpha Polypeptide (CGa). No significant cross-reactivity or interference between Chorionic Gonadotropin Alpha Polypeptide (CGa) and analogues was observed.
Sensitivity:	49.4 pg/mL
Target Details	
Target:	CGA
Alternative Name:	Chorionic Gonadotropin Alpha Polypeptide (CGA Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5652937 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details			
Background:	Gene Name: Chorionic Gonadotropin Alpha Polypeptide Gene Aliases: FSHA, GPHA1, GPHa, HCG, LHA, TSHA, Thyroid Stimulating Hormone Alpha, Follicle-Stimulating Hormone Alpha Subunit, Luteinizing hormone alpha chain, Lutropin alpha chain		
Pathways:	Metabolism of Steroid Hormones and Vitamin D, Thyroid Hormone Synthesis, Hormone Transport, Peptide Hormone Metabolism		
Application Details			
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, ai humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end.		
Assay Time:	2 h		
Plate:	Pre-coated		
Protocol:	This assay employs the competitive inhibition enzyme immunoassay technique. A monoclona antibody specific to Chorionic Gonadotropin Alpha Polypeptide (CGa) has been pre-coated ont a microplate. A competitive inhibition reaction is launched between biotin labeled Chorionic Gonadotropin Alpha Polypeptide (CGa) and unlabeled Chorionic Gonadotropin Alpha Polypeptide (CGa) (Standards or samples) with the pre-coated antibody specific to Chorionic Gonadotropin Alpha Polypeptide (CGa). After incubation the unbound conjugate is washed off. Next, avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. The amount of bound HRP conjugate is reverse proportional to the concentration o Chorionic Gonadotropin Alpha Polypeptide (CGa) in the sample. After addition of the substrate solution, the intensity of color developed is reverse proportional to the concentration of Chorionic Gonadotropin Alpha Polypeptide (CGa) in the sample.		
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Chorionic Gonadotropin Alpha Polypeptide (CGa) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Chorionic Gonadotropin Alpha Polypeptide (CGa) were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100 Intra-Assay: CV<10%		

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5652937 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Application Details	
	Inter-Assay: CV<12%
Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months