

Datasheet for ABIN5653469

CRIM1 ELISA Kit



Overview

Quantity:	96 tests
Target:	CRIM1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1). No significant cross-reactivity or interference between Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) and analogues was observed.
Sensitivity:	0.117 ng/mL

Target Details

Target:	CRIM1
Alternative Name:	Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1 Products)

Target Details Gene Name: Cysteine Rich Transmembrane BMP Regulator 1 Background: Gene Aliases: Cysteine-Rich Motor neuron 1, Chordin-Like, Cysteine-rich repeat-containing protein S52 **Growth Factor Binding** Pathways: **Application Details** Comment: The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. Assay Time: 3 h Plate: Pre-coated Protocol: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1), biotinconjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzymesubstrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Assay Precision: Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level

replicates in each plate. CV(%) = SD/meanX100

Intra-Assay: CV<10%

Inter-Assay: CV<12%

Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) were tested on 3 different plates, 8

Application Details

Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at
	4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant
	pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months