

## Datasheet for ABIN5653606 **DMBT1 ELISA Kit**



Overview

Quantity:	96 tests
Target:	DMBT1
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

## Product Details

Sample Type:	Lung Lavage Fluid, Tissue Homogenate	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Deleted In Malignant Brain Tumors 1 (DMBT1). No significant cross-reactivity or interference between Deleted In Malignant Brain Tumors 1 (DMBT1) and analogues was observed.	
Sensitivity:	0.059 ng/mL	
Target Details		
Target:	DMBT1	

Alternative Name:	Deleted In Malignant Brain Tumors 1 (DMBT1 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5653606 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
Background:	Gene Name: Deleted In Malignant Brain Tumors 1	
	Gene Aliases: GP340, SAG, Muclin, Glycoprotein 340, Hensin, Salivary agglutinin, Surfactant	
	pulmonary-associated D-binding protein	
Pathways:	Activation of Innate immune Response	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Deleted In Malignant Brain	
	Tumors 1 (DMBT1). Standards or samples are then added to the appropriate microtiter plate	
	wells with a biotin-conjugated antibody specific to Deleted In Malignant Brain Tumors 1	
	(DMBT1). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each	
	microplate well and incubated. After TMB substrate solution is added, only those wells that	
	contain Deleted In Malignant Brain Tumors 1 (DMBT1), biotin-conjugated antibody and enzyme-	
	conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by	
	the addition of sulphuric acid solution and the color change is measured	
	spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of Deleted In	
	Malignant Brain Tumors 1 (DMBT1) in the samples is then determined by comparing the O.D. of	
	the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Deleted In Malignant Brain Tumors 1 (DMBT1) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Deleted In Malignant Brain Tumors 1 (DMBT1) were tested on 3 different plates, 8 replicates in	
	each plate. CV(%) = SD/meanX100	
	Intra-Assay: CV<10%	
	Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5653606 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months