

Datasheet for ABIN5653911

Ephrin A5 ELISA Kit



Overview

| Quantity: | 96 tests |
|--------------------------|------------------------|
| Target: | Ephrin A5 (EFNA5) |
| Reactivity: | Human |
| Method Type: | Sandwich ELISA |
| Detection Range: | 0.156 ng/mL - 10 ng/mL |
| Minimum Detection Limit: | 0.156 ng/mL |
| Application: | ELISA |

Product Details

| Sample Type: | Plasma, Serum, Tissue Homogenate |
|--------------------|---|
| Analytical Method: | Quantitative |
| Detection Method: | Colorimetric |
| Specificity: | This assay has high sensitivity and excellent specificity for detection of Ephrin A5 (EFNA5). No significant cross-reactivity or interference between Ephrin A5 (EFNA5) and analogues was observed. |
| Sensitivity: | 0.055 ng/mL |

Target Details

| Target: | Ephrin A5 (EFNA5) |
|-------------------|----------------------------|
| Alternative Name: | Ephrin A5 (EFNA5 Products) |

Target Details

| Background: | Gene Name: Ephrin A5 |
|---------------------|---|
| | Gene Aliases: AF1, EFL5, EPLG7, GLC1M, LERK7, RAGS, AL-1, EPH-related receptor tyrosine |
| | kinase ligand 7 |
| Gene ID: | 1946 |
| UniProt: | P52803 |
| Pathways: | RTK Signaling |
| Application Details | |
| Comment: | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than |
| | 5 % within the expiration date under appropriate storage condition. To minimize extra influence |
| | on the performance, operation procedures and lab conditions, especially room temperature, air |
| | humidity, incubator temperature should be strictly controlled. It is also strongly suggested that |
| | the whole assay is performed by the same operator from the beginning to the end. |
| Assay Time: | 3 h |
| Plate: | Pre-coated |
| Protocol: | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate |
| | provided in this kit has been pre-coated with an antibody specific to Ephrin A5 (EFNA5). |
| | Standards or samples are then added to the appropriate microtiter plate wells with a biotin- |
| | conjugated antibody specific to Ephrin A5 (EFNA5). Next, Avidin conjugated to Horseradish |
| | Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution |
| | is added, only those wells that contain Ephrin A5 (EFNA5), biotin-conjugated antibody and |
| | enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is |
| | terminated by the addition of sulphuric acid solution and the color change is measured |
| | spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Ephrin A5 |
| | (EFNA5) in the samples is then determined by comparing the O.D. of the samples to the |
| | standard curve. |
| Assay Precision: | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level |
| | Ephrin A5 (EFNA5) were tested 20 times on one plate, respectively |
| | Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level |
| | Ephrin A5 (EFNA5) were tested on 3 different plates, 8 replicates in each plate. CV(%) = |
| | |
| | SD/meanX100 |
| | SD/meanX100 Intra-Assay: CV<10% |

Application Details

| Restrictions: | For Research Use only |
|------------------|--|
| Handling | |
| Handling Advice: | The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all |
| | steps should be completed without interruption. |
| Storage: | 4 °C,-20 °C |
| Storage Comment: | -20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles. |
| Expiry Date: | 4-8 months |