

Datasheet for ABIN5654047 FAM3A ELISA Kit



Overview

| Quantity: | 96 tests |
|--------------------------|--------------------------|
| Target: | FAM3A |
| Reactivity: | Human |
| Method Type: | Sandwich ELISA |
| Detection Range: | 15.62 pg/mL - 1000 pg/mL |
| Minimum Detection Limit: | 15.62 pg/mL |
| Application: | ELISA |

Product Details

| Sample Type: | Cell Lysate, Tissue Homogenate | |
|--------------------|---|--|
| Analytical Method: | Quantitative | |
| Detection Method: | Colorimetric | |
| Specificity: | This assay has high sensitivity and excellent specificity for detection of Family With Sequence Similarity 3, Member A (FAM3A). No significant cross-reactivity or interference between Family With Sequence Similarity 3, Member A (FAM3A) and analogues was observed. | |
| Sensitivity: | 26.1 pg/mL | |
| Target Details | | |
| Target: | FAM3A | |

| Alternative Name: | Family With Sequence Similarity 3, Member A (FAM3A Products) |
|-------------------|--|
|-------------------|--|

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5654047 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details

Background:

Gene Name: Family With Sequence Similarity 3, Member A Gene Aliases: XAP-7, Cytokine-like protein 2-19

Application Details

| Comment: | The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than |
|------------------|---|
| | 5 % within the expiration date under appropriate storage condition. To minimize extra influence |
| | on the performance, operation procedures and lab conditions, especially room temperature, air |
| | humidity, incubator temperature should be strictly controlled. It is also strongly suggested that |
| | the whole assay is performed by the same operator from the beginning to the end. |
| Assay Time: | 3 h |
| Plate: | Pre-coated |
| Protocol: | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate |
| | provided in this kit has been pre-coated with an antibody specific to Family With Sequence |
| | Similarity 3, Member A (FAM3A). Standards or samples are then added to the appropriate |
| | microtiter plate wells with a biotin-conjugated antibody specific to Family With Sequence |
| | Similarity 3, Member A (FAM3A). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is |
| | added to each microplate well and incubated. After TMB substrate solution is added, only those |
| | wells that contain Family With Sequence Similarity 3, Member A (FAM3A), biotin-conjugated |
| | antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate |
| | reaction is terminated by the addition of sulphuric acid solution and the color change is |
| | measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of |
| | Family With Sequence Similarity 3, Member A (FAM3A) in the samples is then determined by |
| | comparing the O.D. of the samples to the standard curve. |
| Assay Precision: | Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level |
| | Family With Sequence Similarity 3, Member A (FAM3A) were tested 20 times on one plate, |
| | respectively |
| | Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level |
| | Family With Sequence Similarity 3, Member A (FAM3A) were tested on 3 different plates, 8 |
| | replicates in each plate. CV(%) = SD/meanX100 |
| | Intra-Assay: CV<10% |
| | Inter-Assay: CV<12% |
| Restrictions: | For Research Use only |
| | |

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5654047 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

| Н | land | lina |
|----|------|------|
| 11 | ana | mig |

| Handling Advice: | The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption. |
|------------------|--|
| Storage: | 4 °C,-20 °C |
| Storage Comment: | -20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles. |
| Expiry Date: | 4-8 months |