

Datasheet for ABIN5654196 FGF11 ELISA Kit



Overview

Quantity:	96 tests
Target:	FGF11
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	15.62 pg/mL - 1000 pg/mL
Minimum Detection Limit:	15.62 pg/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Fibroblast Growth Factor 11 (FGF11). No significant cross-reactivity or interference between Fibroblast Growth Factor 11 (FGF11) and analogues was observed.
Sensitivity:	7.2 pg/mL
Target Details	
Target:	FGF11
Alternative Name:	Fibroblast Growth Factor 11 (FGF11 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5654196 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target	Details
--------	---------

Background:

Gene Name: Fibroblast Growth Factor 11 Gene Aliases: FHF3, Fibroblast Growth Factor Homologous Factor 3

Application Details

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
Handling	
Restrictions:	For Research Use only
	Inter-Assay: CV<12%
	Intra-Assay: CV<10%
	plate. CV(%) = SD/meanX100
	Fibroblast Growth Factor 11 (FGF11) were tested on 3 different plates, 8 replicates in each
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Fibroblast Growth Factor 11 (FGF11) were tested 20 times on one plate, respectively
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	by comparing the O.D. of the samples to the standard curve.
	The concentration of Fibroblast Growth Factor 11 (FGF11) in the samples is then determined
	and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm.
	color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution
	11 (FGF11), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in
	After TMB substrate solution is added, only those wells that contain Fibroblast Growth Factor
	conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated.
	a biotin-conjugated antibody specific to Fibroblast Growth Factor 11 (FGF11). Next, Avidin
	11 (FGF11). Standards or samples are then added to the appropriate microtiter plate wells wit
	provided in this kit has been pre-coated with an antibody specific to Fibroblast Growth Factor
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
Plate:	Pre-coated
Assay Time:	3 h
	the whole assay is performed by the same operator from the beginning to the end.
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	on the performance, operation procedures and lab conditions, especially room temperature, ai
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5654196 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

	specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months