

Datasheet for ABIN5654342 FCN1 ELISA Kit



Overview

Quantity:	96 tests
Target:	FCN1
Reactivity:	Mouse
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

Product Details

Alternative Name:

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Ficolin 1 (FCN1). No significant cross-reactivity or interference between Ficolin 1 (FCN1) and analogues was observed.
Sensitivity:	0.056 ng/mL
Target Details	
Target:	FCN1

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5654342 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Ficolin 1 (FCN1 Products)

Background:	Gene Name: Ficolin 1
	Gene Aliases: FCNM, Collagen/fibrinogen domain-containing protein 1, Ficolin-alpha, Ficolin-A,
	M-ficolin
Gene ID:	14133
UniProt:	070165
Pathways:	Complement System
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to Ficolin 1 (FCN1). Standards
	or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated
	antibody specific to Ficolin 1 (FCN1). Next, Avidin conjugated to Horseradish Peroxidase (HRP)
	is added to each microplate well and incubated. After TMB substrate solution is added, only
	those wells that contain Ficolin 1 (FCN1), biotin-conjugated antibody and enzyme-conjugated
	Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the
	addition of sulphuric acid solution and the color change is measured spectrophotometrically at
	a wavelength of 450nm \pm 10nm. The concentration of Ficolin 1 (FCN1) in the samples is then
	determined by comparing the O.D. of the samples to the standard curve.

Restrictions:

Assay Precision:

For Research Use only

Intra-Assay: CV<10% Inter-Assay: CV<12%

SD/meanX100

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5654342 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Ficolin 1 (FCN1) were tested 20 times on one plate, respectively

Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level

Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level

Ficolin 1 (FCN1) were tested on 3 different plates, 8 replicates in each plate. CV(%) =

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months