

Datasheet for ABIN5654699 **GPX3 ELISA Kit**



Overview

Quantity:	96 tests
Target:	GPX3
Reactivity:	Rat
Method Type:	Sandwich ELISA
Detection Range:	0.781 ng/mL - 50 ng/mL
Minimum Detection Limit:	0.781 ng/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum	
Analytical Method:	Quantitative	
Detection Method:	Colorimetric	
Specificity:	This assay has high sensitivity and excellent specificity for detection of Glutathione Peroxidase 3, Plasma (GPX3). No significant cross-reactivity or interference between Glutathione Peroxidase 3, Plasma (GPX3) and analogues was observed.	
Sensitivity:	0.37 ng/mL	
Target Details		
Target:	GPX3	
Alternative Name:	Glutathione Peroxidase 3, Plasma (GPX3 Products)	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5654699 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details	
Background:	Gene Name: Glutathione Peroxidase 3, Plasma
	Gene Aliases: GSHPX3, GPx-P, GSHPx-P, GPXP, Plasma glutathione peroxidase, Extracellular glutathione peroxidase
Pathways:	Thyroid Hormone Synthesis
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 $\%$ within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to Glutathione Peroxidase 3,
	Plasma (GPX3). Standards or samples are then added to the appropriate microtiter plate wells
	with a biotin-conjugated antibody specific to Glutathione Peroxidase 3, Plasma (GPX3). Next,
	Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and
	incubated. After TMB substrate solution is added, only those wells that contain Glutathione
	Peroxidase 3, Plasma (GPX3), biotin-conjugated antibody and enzyme-conjugated Avidin will
	exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of
	sulphuric acid solution and the color change is measured spectrophotometrically at a
	wavelength of 450nm \pm 10nm. The concentration of Glutathione Peroxidase 3, Plasma (GPX3)
	in the samples is then determined by comparing the O.D. of the samples to the standard curve
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	Glutathione Peroxidase 3, Plasma (GPX3) were tested 20 times on one plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Glutathione Peroxidase 3, Plasma (GPX3) were tested on 3 different plates, 8 replicates in each
	plate. CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%
Restrictions:	For Research Use only

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5654699 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months