

# Datasheet for ABIN5655117

### **HEXB ELISA Kit**



#### Overview

Quantity:	96 tests
Target:	HEXB
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.156 ng/mL - 10 ng/mL
Minimum Detection Limit:	0.156 ng/mL
Application:	ELISA

### **Product Details**

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Hexosaminidase B Beta (HEXb). No significant cross-reactivity or interference between Hexosaminidase B Beta (HEXb) and analogues was observed.
Sensitivity:	0.058 ng/mL

# **Target Details**

Target:	HEXB
Alternative Name:	Hexosaminidase B Beta (HEXB Products)

# Target Details

Background:	Gene Name: Hexosaminidase B Beta
	Gene Aliases: Hex-B, ENC-1AS, HCC-7, Beta-N-acetylhexosaminidase subunit beta, Cervical
	cancer proto-oncogene 7 protein, N-acetyl-beta-glucosaminidase subunit beta
Gene ID:	3074
UniProt:	P07686
Pathways:	Sensory Perception of Sound, Glycosaminoglycan Metabolic Process
Application Details	
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
	5 % within the expiration date under appropriate storage condition. To minimize extra influence
	on the performance, operation procedures and lab conditions, especially room temperature, air
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
	the whole assay is performed by the same operator from the beginning to the end.
Assay Time:	3 h
Plate:	Pre-coated
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
	provided in this kit has been pre-coated with an antibody specific to Hexosaminidase B Beta
	(HEXb). Standards or samples are then added to the appropriate microtiter plate wells with a
	biotin-conjugated antibody specific to Hexosaminidase B Beta (HEXb). Next, Avidin conjugated
	to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB
	substrate solution is added, only those wells that contain Hexosaminidase B Beta (HEXb),
	biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The
	enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the colo
	change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The
	concentration of Hexosaminidase B Beta (HEXb) in the samples is then determined by
	comparing the O.D. of the samples to the standard curve.
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
	Hexosaminidase B Beta (HEXb) were tested 20 times on one plate, respectively
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
	Hexosaminidase B Beta (HEXb) were tested on 3 different plates, 8 replicates in each plate.
	CV(%) = SD/meanX100
	Intra-Assay: CV<10%
	Inter-Assay: CV<12%

# **Application Details**

Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at
	4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant
	pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months