

## Datasheet for ABIN5655369

### **INPP5D ELISA Kit**



#### Overview

Quantity:	96 tests
Target:	INPP5D
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	31.25 pg/mL - 2000 pg/mL
Minimum Detection Limit:	31.25 pg/mL
Application:	ELISA

#### **Product Details**

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D). No significant cross-reactivity or interference between Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D) and analogues was observed.
Sensitivity:	13.5 pg/mL

## **Target Details**

Target:	INPP5D
Alternative Name:	Inositol Polyphosphate-5-Phosphatase 145kDa (INPP5D Products)

# Target Details Gene Name: Inositol Polyphosphate-5-Phosphatase 145 kDa Background: Gene Aliases: SHIP, SHIP1, SIP-145, hp51CN, Phosphatidylinositol-3,4,5-Trisphosphate 5-Phosphatase 1, SH2 domain-containing inositol 5'-phosphatase 1 TCR Signaling, BCR Signaling, Warburg Effect Pathways: **Application Details** Comment: The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than 5 % within the expiration date under appropriate storage condition. To minimize extra influence on the performance, operation procedures and lab conditions, especially room temperature, air humidity, incubator temperature should be strictly controlled. It is also strongly suggested that the whole assay is performed by the same operator from the beginning to the end. Assay Time: 3 h Plate: Pre-coated Protocol: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D). Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody specific to Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D). Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Inositol Polyphosphate-5-Phosphatase 145 kDa (INPP5D) in the samples is then determined by comparing the O.D. of the samples to the standard curve. Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level Assay Precision: Inositol Polyphosphate-5-Phosphatase 145kDa (INPP5D) were tested 20 times on one plate, respectively Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level Inositol Polyphosphate-5-Phosphatase 145kDa (INPP5D) were tested on 3 different plates, 8

replicates in each plate. CV(%) = SD/meanX100

Intra-Assay: CV<10%

Inter-Assay: CV<12%

## **Application Details**

Restrictions:	For Research Use only
Handling	
Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and
	specimen samples should be assayed in duplicate. Once the procedure has been started, all
	steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at
	4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant
	pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months