

Datasheet for ABIN5656124 IARS ELISA Kit



Overview

Quantity:	96 tests
Target:	IARS
Reactivity:	Human
Method Type:	Sandwich ELISA
Detection Range:	0.312 ng/mL - 20 ng/mL
Minimum Detection Limit:	0.312 ng/mL
Application:	ELISA

Product Details

Sample Type:	Cell Lysate, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Isoleucyl tRNA Synthetase (IARS). No significant cross-reactivity or interference between Isoleucyl tRNA Synthetase (IARS) and analogues was observed.
Sensitivity:	0.114 ng/mL
Target Details	
Target:	IARS

Alternative Name:	Isoleucyl tRNA Synthetase (IARS Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5656124 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Target Details		
Background:	Gene Name: Isoleucyl tRNA Synthetase	
	Gene Aliases: IARS1, ILRS, OJ, IleRS, IRS, Isoleucine tRNA Ligase 1,Cytoplasmic	
Gene ID:	3376	
UniProt:	P41252	
Application Details		
Comment:	The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than	
	5 % within the expiration date under appropriate storage condition. To minimize extra influence	
	on the performance, operation procedures and lab conditions, especially room temperature, air	
	humidity, incubator temperature should be strictly controlled. It is also strongly suggested that	
	the whole assay is performed by the same operator from the beginning to the end.	
Assay Time:	3 h	
Plate:	Pre-coated	
Protocol:	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate	
	provided in this kit has been pre-coated with an antibody specific to Isoleucyl tRNA Synthetase	
	(IARS). Standards or samples are then added to the appropriate microtiter plate wells with a	
	biotin-conjugated antibody specific to Isoleucyl tRNA Synthetase (IARS). Next, Avidin	
	conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated.	
	After TMB substrate solution is added, only those wells that contain Isoleucyl tRNA Synthetase	
	(IARS), biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color	
	The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the	
	color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The	
	concentration of Isoleucyl tRNA Synthetase (IARS) in the samples is then determined by	
	comparing the O.D. of the samples to the standard curve.	
Assay Precision:	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level	
	Isoleucyl tRNA Synthetase (IARS) were tested 20 times on one plate, respectively	
	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level	
	Isoleucyl tRNA Synthetase (IARS) were tested on 3 different plates, 8 replicates in each plate.	
	CV(%) = SD/meanX100	
	Intra-Assay: CV<10%	
	Inter-Assay: CV<12%	
Restrictions:	For Research Use only	

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5656124 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months