

Datasheet for ABIN5657303 NPX1 ELISA Kit



Overview

Quantity:	96 tests
Target:	NPX1
Reactivity:	Rat
Method Type:	Sandwich ELISA
Detection Range:	15.62 pg/mL - 1000 pg/mL
Minimum Detection Limit:	15.62 pg/mL
Application:	ELISA

Product Details

Sample Type:	Plasma, Serum, Tissue Homogenate
Analytical Method:	Quantitative
Detection Method:	Colorimetric
Specificity:	This assay has high sensitivity and excellent specificity for detection of Neuronal Pentraxin I (NPTX1). No significant cross-reactivity or interference between Neuronal Pentraxin I (NPTX1) and analogues was observed.
Sensitivity:	5.5 pg/mL
Target Details	
Target:	NPX1
Alternative Name:	Neuronal Pentraxin I (NPX1 Products)

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/3 | Product datasheet for ABIN5657303 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Gene Name: Neuronal Pentraxin I
Gene Aliases: NARP1, NP-I, NP1
Carbohydrate Homeostasis, Regulation of Cell Size, Signaling Events mediated by VEGFR1 and
VEGFR2, Smooth Muscle Cell Migration, Platelet-derived growth Factor Receptor Signaling,
VEGFR1 Specific Signals, SARS-CoV-2 Protein Interactome
The stability of kit is determined by the loss rate of activity. The loss rate of this kit is less than
5 % within the expiration date under appropriate storage condition. To minimize extra influence
on the performance, operation procedures and lab conditions, especially room temperature, air
humidity, incubator temperature should be strictly controlled. It is also strongly suggested that
the whole assay is performed by the same operator from the beginning to the end.
3 h
Pre-coated
The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate
provided in this kit has been pre-coated with an antibody specific to Neuronal Pentraxin I
(NPTX1). Standards or samples are then added to the appropriate microtiter plate wells with a
biotin-conjugated antibody specific to Neuronal Pentraxin I (NPTX1). Next, Avidin conjugated to
Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB
substrate solution is added, only those wells that contain Neuronal Pentraxin I (NPTX1), biotin-
conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-
substrate reaction is terminated by the addition of sulphuric acid solution and the color change
is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of
Neuronal Pentraxin I (NPTX1) in the samples is then determined by comparing the O.D. of the
samples to the standard curve.
Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level
Neuronal Pentraxin I (NPTX1) were tested 20 times on one plate, respectively
Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level
Neuronal Pentraxin I (NPTX1) were tested on 3 different plates, 8 replicates in each plate. CV(%)
= SD/meanX100
Intra-Assay: CV<10%
Inter-Assay: CV<12%
For Research Use only

Order at www.antibodies-online.com | www.antikoerper-online.de | www.anticorps-enligne.fr | www.antibodies-online.cn International: +49 (0)241 95 163 153 | USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 2/3 | Product datasheet for ABIN5657303 | 07/25/2024 | Copyright antibodies-online. All rights reserved.

Н	land	lina
11	ana	mig

Handling Advice:	The Stop Solution is acidic. Do not allow to contact skin or eyes. Calibrators, controls and specimen samples should be assayed in duplicate. Once the procedure has been started, all steps should be completed without interruption.
Storage:	4 °C,-20 °C
Storage Comment:	-20°C. Bring all reagents to room temperature before beginning test. The kit may be stored at 4°C for immediate use within two days upon arrival. Reseal any unused strips with desiccant pack. Minimize freeze/thaw cycles.
Expiry Date:	4-8 months